


Microcentrifuge Tube, Non-sterile 마이크로 튜브

Cat. No.	Model	규격	Color	PK Price ₩	CS Price ₩
T16-510-010	MT61015	1.5 ml	Clear	500ea 11,000 10pk	99,000

- 재질 : PP, -80 ~ 121°C 내열, Autoclavable
- Max RCF : 16000×g
- DNase, RNase & Pyrogen Free
- LK Labkorea



Boil-Proof Microtube 마이크로 튜브, E-P Tube

Cat. No.	Model	용량	Cap	PK Price ₩
T16-46-601	T50000-N	0.5 ml	Clear	1000ea 24,600
T16-46-619	T50010-N	1.5 ml		1000ea 30,500
T16-46-637	T50020-N	2.0 ml		500ea 13,000

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- Max RCF : 20000×g
- DNase, RNase & Pyrogen Free
- Tarsons



Microcentrifuge Tube / Microtube, PCR-suitable 마이크로 튜브, BIO-CERT® PCR QUALITY

- 재질 : High-purity PP, Max. 135°C 내열, Autoclavable
- RNase, DNA & Pyrogen Free (acc. to LAL-test)
- PCR Inhibitor Free
- High Cleanliness : PCR-setup, Storage of PCR Sample 등
- Brand GmBH

with Easy-to-open Snap-on Lid


Cat. No.	Model	용량	색상	Max. RCF	PK Price ₩	
T16-59-421	780507	0.5 ml ①	Natural	10000	1000ea 85,500	
T16-59-471	780420	1.5 ml ②			500ea 36,100	
T16-59-473	780421		Yellow	500ea 49,300		
T16-59-475	780422		Blue	500ea 49,300		
T16-59-477	780423		Green	500ea 49,300		
T16-59-479	780424		Orange	500ea 49,300		
T16-59-481	780425		Amber	500ea 49,300		
T16-59-461	780550		2.0 ml ③	Natural	25000	500ea 56,600
T16-59-467	780555		5.0 ml ④			250ea 100,700

- 권장 Storage Temp. : up to -20°C
- Amber Color : Suitable for Light-Sensitive Reagent

with Easy-to-open Snap-on Lid & Lid Locking

Cat. No.	Model	용량	색상	Max. RCF	PK Price ₩
T16-59-490	780536	0.5 ml ⑤	Natural	30000	500ea 45,400
T16-59-492	780540	1.5 ml ⑥			1000ea 78,800
T16-59-494	780546	2.0 ml ⑦			500ea 47,600

- 권장 Storage Temp. : up to -80°C
- Sample 가열 시에 뚜껑이 열리는 것을 방지하는 이중 Lock 구조



Microtube, without Lid 마이크로 튜브

Cat. No.	Model	용량	색상	Max. RCF	PK Price ₩
T16-59-451	780505	1.5 ml	Transparent	6000	2000ea 111,400

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- Size : Ø11 × h39.5 mm
- Non-Sterile
- Brand GmBH

