

### Microtube with Screw Cap & Silicon Seal, Sterile 멸균 마이크로 튜브, BIO-CERT® STERILE

- 재질 : PP Tube, PP Cover Thread, Silicone Seal
- Temp. Range : -196 °C ~ 121 °C, Autoclavable
- Max. RCF : 17000 g (with base, Flat) / 20000 g (without base, Conical)
- Non-mutagenic, Non-cytotoxic, DNA, DNase, RNase & Endotoxins Free
- Cap Insert는 1134 Page 참조
- Brand GmBH




#### without Lid Closure, to Protect against Evaporation ①

Cat. No.	Model	용량	눈금	Bottom	PK Price ₩
C25-59-201	780750	0.5 ml	No	Flat	500ea 367,200
C25-59-205	780751				500ea 367,200
C25-59-209	780752	1.5 ml	Yes	Conical	500ea 367,200
C25-59-213	780753			Flat	500ea 367,200
C25-59-217	780754	2.0 ml		Conical	500ea 367,200
				Flat	500ea 367,200

#### with Lid Closure, for Storage in Gaseous of LN<sub>2</sub> ②

Cat. No.	Model	용량	눈금	Bottom	PK Price ₩
C25-59-221	780755	0.5 ml	No	Flat	500ea 502,800
C25-59-225	780756				500ea 502,800
C25-59-229	780757	1.5 ml	Yes	Conical	500ea 502,800
C25-59-223	780758			Flat	500ea 502,800
C25-59-227	780759	2.0 ml		Conical	500ea 502,800
				Flat	500ea 502,800

- Temper-evident screw cap guarantees Uncontamination Sample



### Cryotube / Cryogenic Vial, External Thread 멸균 냉동 튜브 / 멸균 냉동 바이알, Sterile

Cat. No.	Model	용량	눈금	Bottom	CS Price ₩
T16-59-231	114830	1.2 ml	0.1 ml	Self-standing	1000ea 684,500
T16-59-233	114831	2.0 ml		Round	1000ea 684,500
T16-59-235	114832			1000ea 684,500	
T16-59-237	114833	3.0 ml		Self-standing	1000ea 732,400
T16-59-239	114834	4.0 ml		1000ea 793,300	
T16-59-241	114835	5.0 ml		1000ea 837,600	

- 재질 : PP Tube & Cover, Silicon Seal, -196 °C ~ 121 °C, Autoclavable
- -196 °C까지의 Long-term Storage에 맞맞음
- DNA, DNase, RNase & Endotoxins Free
- Brand GmBH



### Cryotube / Cryogenic Vial, Internal Thread 멸균 냉동 튜브 / 멸균 냉동 바이알, Sterile

Cat. No.	Model	용량	눈금	Bottom	CS Price ₩
T16-59-245	114840	1.2 ml	0.1 ml	Self-standing	1000ea 684,500
T16-59-247	114841	2.0 ml		Round	1000ea 684,500
T16-59-249	114842			1000ea 684,500	
T16-59-251	114843	4.0 ml		Self-standing	1000ea 911,600
T16-59-253	114844			1000ea 1,038,300	
T16-59-255	114845	5.0 ml		Round	1000ea 1,040,800

- 재질 : PP Tube & Cover, Silicon Seal, -196 °C ~ 121 °C, Autoclavable
- -196 °C까지의 Long-term Storage에 맞맞음
- DNA, DNase, RNase & Endotoxins Free
- Brand GmBH



#### Cap Insert for "Brand" Cryogenic Vial

Cat. No.	Model	Color	PK Price ₩
T16-59-901	114850	White	500ea 34,700
T16-59-903	114851	Blue	500ea 34,700
T16-59-905	114852	Red	500ea 34,700
T16-59-907	114853	Green	500ea 34,700
T16-59-909	114854	Yellow	500ea 34,700

- 재질 : PP • Fit for all size

