


Culture Tube with Screw Cap 플라스틱 캡 시험관, 캡 컬처 튜브

Cat. No.	Model	용량	Size	PK Price ₩	CS Price ₩
T16-660-131	PCT1610	PP Tube	∅16×h100 mm	100ea 14,200	10pk 118,900

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- Natural Color, Assembled Screw Cap 포함
- LK Labkorea



Sterile Culture Tube with Dual-position Cap 듀얼 캡 컬처 튜브, 멸균



- 튜브 재질 : PP, -10 ~ 135°C 내열, Autoclavable,
- 튜브 재질 : PS, -10 ~ 80°C 내열, Non-autoclavable
- Dual-position Ribbed PE Cap, Aerobic 또는 Anaerobic Culture용
- LK Labkorea

PP Test Tube, Sterile, 반투명

Cat. No.	Model	용량	눈금	Size ∅×h	PK Price ₩	CS Price ₩
T16-660-151	CT8501	5 ml	No	12×75 mm	25ea 2,900	40pk 97,800
T16-660-155	CT8505	15 ml		17×100 mm	25ea 5,300	40pk 183,000
T16-660-161	CT8521	5 ml	Yes	12×75 mm	25ea 3,600	40pk 123,100
T16-660-165	CT8525	15 ml		17×100 mm	25ea 6,000	40pk 208,300

PS Test Tube, Sterile, 투명

Cat. No.	Model	용량	눈금	Size ∅×h	PK Price ₩	CS Price ₩
T16-660-171	CT8541	5 ml	No	12×75 mm	25ea 2,900	40pk 99,600
T16-660-175	CT8545	15 ml		17×100 mm	25ea 5,300	40pk 183,000
T16-660-181	CT8561	5 ml	Yes	12×75 mm	25ea 3,400	40pk 117,700
T16-660-185	CT8565	15 ml		17×100 mm	25ea 6,000	40pk 208,300



PFA Test Tube with Screw Cap PFA 캡 시험관

Cat. No.	Model	용량	Size	EA Price ₩
T16-64-781	F-1181-002	15 ml	∅16×h110 mm	263,700
T16-64-785	F-1181-003	26 ml	∅22×h112 mm	347,300
T16-64-789	F-1181-004	32 ml	∅20×h150 mm	207,300
T16-64-793	F-1181-005	60 ml	∅26×h164 mm	418,300

- 재질 : PFA (Perfluoroalkoxy)
- -196 ~ +250°C 내열
- Flon Industry



PTFE Centrifuge Test Tube PTFE 시험관 / 원심관

- 재질 : PTFE (Polytetrafluoroethylene)
- -267 ~ +280°C 내열
- Cowie

Round Bottom (입구 : Lip)

Cat. No.	Model	용량	Size	EA Price ₩
T16-450-211	012.108	8 ml	∅12×100 mm	19,000
T16-450-215	012.113	13 ml	∅16×100 mm	18,900
T16-450-219	012.118	17 ml	∅18×100 mm	25,900
T16-450-223	012.150	50 ml	∅25×145 mm	33,100
T16-450-227	012.180	85 ml	∅40×100 mm	33,100

Round Bottom with Screw Cap

Cat. No.	Model	용량	Size	EA Price ₩
T16-450-111	012.1108	8 ml	∅12×100 mm	27,100
T16-450-115	012.1118	17 ml	∅18×100 mm	32,500
T16-450-119	012.1150	50 ml	∅25×140 mm	55,900



Conical Bottom (입구 : Lip)

Cat. No.	Model	용량	Size	EA Price ₩
T16-450-411	012.213	13 ml	∅16×110 mm	18,600
T16-450-415	012.215	18 ml	∅18×120 mm	24,100
T16-450-419	012.245	45 ml	∅25×150 mm	40,100

Conical Bottom with Screw Cap

Cat. No.	Model	용량	Size	EA Price ₩
T16-450-311	012.1213	13 ml	∅16×110 mm	28,400
T16-450-315	012.1215	18 ml	∅18×120 mm	35,200
T16-450-319	012.1245	45 ml	∅25×150 mm	54,000