


I Tungsten Boat 텡스텐 보트

Cat. No.	Model	높이 (h)	폭 (w)	두께 (t)	EA Price ₩
B07-142-100	LTB201	2.0 mm	8 mm	0.2 t	5,100
B07-142-110	LTB202	2.2 mm	10 mm		6,400
B07-142-120	LTB203	3.0 mm	15 mm		10,200
B07-142-130	LTB204	3.5 mm	20 mm		10,200
B07-142-140	LTB205	4.0 mm	25 mm		11,900
B07-142-150	LTB301	2.0 mm	8 mm	0.3 t	6,800
B07-142-160	LTB302	2.2 mm	10 mm		6,800
B07-142-170	LTB303	3.0 mm	15 mm		10,200
B07-142-180	LTB304	4.0 mm	20 mm		13,500
B07-142-190	LTB305	4.5 mm	25 mm		16,900
B07-142-200	LTB306	5.0 mm	30 mm	18,600	
B07-142-210	LTB401	1.0 mm	8 mm	0.5 t	11,000
B07-142-220	LTB402	2.2 mm	10 mm		11,000
B07-142-230	LTB403	3.5 mm	15 mm		16,100
B07-142-240	LTB404	4.0 mm	20 mm		22,000
B07-142-250	LTB405	5.0 mm	25 mm		25,400
B07-142-260	LTB406	5.3 mm	30 mm		28,700

- 전체 길이는 100 mm, 홈 길이는 50 mm
- LK Labkorea



I Porcelain Weighing Boat 자제 평량 보트

Cat. No.	Model	Size	EA Price ₩	PK Price ₩
B07-128-115	252/0	w37 × d17 mm	3,600	10ea 34,900
B07-128-125	252/1	w53 × d23 mm	4,000	10ea 37,300
B07-128-135	252/2	w64 × d29 mm	4,100	10ea 38,400

- 재질 : Porcelain, 1,000°C 내열
- Jipo / LabCermic®



- Paper Weighing Boats는 59 Page 참조
- Plastic Weighing Boats는 61 Page 참조

