

Screwcap Culture Tube PP 캡 시험관

Cat. No.	Model	용량	Size (Ø×h)	Cap	PK Price ₩	
T16-59-211	113931	7 ml	12×100 mm	GL14	100ea	161,900
T16-59-213	113935	15 ml	16×100 mm	GL18	100ea	144,200
T16-59-215	113941	23 ml	16×160 mm		100ea	206,100
T16-59-217	113943	36 ml	18×180 mm		100ea	253,900

- Soda-lime Glass, Autoclavable
- 두께 1t의 보급형 캡 시험관
- PP Screwcap with White Elastomer Seal
- Brand GmBH



Screwcap Tube 캡 시험관, LukeGL®

- Borosilicate Glass 3.3
- LukeGL®

with Rubber Lined Screwcap, Culture 용, 멸균가능(Autoclavable)

Cat. No.	Model	용량	Size (OD×h)	Cap	PK Price ₩		CS Price ₩	
T16-256-621	TU-RC11	11 ml	16×100 mm	15-425	10ea	26,600	5pk	121,200
T16-256-628	TU-RC14	14 ml	16×125 mm		10ea	27,500	5pk	125,200
T16-256-635	TU-RC18	18 ml	16×150 mm		10ea	29,300	5pk	133,000
T16-256-642	TU-RC25	25 ml	18×150 mm	18-400	10ea	33,000	5pk	150,200
T16-256-649	TU-RC30	30 ml	18×180 mm		10ea	34,700	5pk	158,100
T16-256-656	TU-RC32	32 ml	20×150 mm		10ea	36,500	5pk	165,900
T16-256-663	TU-RC40	40 ml	20×180 mm	24-400	10ea	38,200	5pk	173,700
T16-256-670	TU-RC50	50 ml	26×150 mm		5ea	27,500	5pk	124,800
T16-256-677	TU-RC72	72 ml	26×200 mm		5ea	30,000	5pk	136,500

- 멸균시 반드시 뚜껑을 풀어 내압이 빠져갈 수 있도록 해야함

with Teflon Lined Cap, Solvent 용 (Non-Autoclavable)

Cat. No.	Model	용량	Size (OD×h)	Cap	PK Price ₩		CS Price ₩	
T16-256-721	TU-TC11	11 ml	16×100 mm	15-425	10ea	25,300	5pk	115,200
T16-256-727	TU-TC14	14 ml	16×125 mm		10ea	29,400	5pk	132,500
T16-256-733	TU-TC18	18 ml	16×150 mm		10ea	29,800	5pk	133,600
T16-256-739	TU-TC25	25 ml	18×150 mm	18-400	10ea	36,500	5pk	164,300
T16-256-745	TU-TC30	30 ml	18×180 mm		10ea	37,300	5pk	168,100
T16-256-751	TU-TC32	32 ml	20×150 mm		10ea	37,700	5pk	170,100
T16-256-757	TU-TC40	40 ml	20×180 mm	24-400	10ea	38,600	5pk	173,700
T16-256-763	TU-TC50	50 ml	26×150 mm		5ea	28,200	5pk	128,300
T16-256-769	TU-TC72	72 ml	26×200 mm		5ea	30,800	5pk	140,100

- 캡 내부가 Teflon으로 되어진 것으로 최상의 내화학성을 가짐



with Teflon Lined Cap

- Optional Tube Cap은 1249 ~ 1251 Page 참조



PBT Screwcap Tube PBT 캡 시험관, LukeGL®

Cat. No.	Model	용량	Size (od×h)	Cap	PK Price ₩		CS Price ₩	
T16-256-322	TV-RC11	11 ml	16×100	GL14	5ea	32,200	5pk	146,500
T16-256-326	TV-RC14	14 ml	16×125		5ea	32,700	5pk	148,600
T16-256-330	TV-RC18	18 ml	16×150		5ea	33,500	5pk	152,400
T16-256-334	TV-RC25	25 ml	18×150	GL18	5ea	38,600	5pk	175,200
T16-256-338	TV-RC30	30 ml	18×180		5ea	39,400	5pk	179,200
T16-256-342	TV-RC32	32 ml	20×150		5ea	40,200	5pk	183,000
T16-256-346	TV-RC40	40 ml	20×180	GL25	5ea	41,100	5pk	187,000
T16-256-350	TV-RC50	50 ml	26×150		5ea	56,100	2pk	102,300
T16-256-354	TV-RC72	72 ml	26×200		5ea	58,800	2pk	107,000
T16-256-358	TV-RC100	100 ml	30×200	GL32	5ea	67,300	2pk	122,700
T16-256-362	TV-RC130	130 ml	30×250		5ea	70,700	2pk	127,300
T16-256-366	TV-RC101	100 ml	38×150		1ea	25,300	5pk	113,900
T16-256-370	TV-RC150	150 ml	38×200	GL32	1ea	26,400	5pk	118,600
T16-256-374	TV-RC200	200 ml	38×250		1ea	30,500	5pk	137,300
T16-256-378	TV-RC250	250 ml	38×300		1ea	32,800	5pk	147,600

- Borosilicate Glass 3.3
- Autoclavable
- PTFE Coated Silicon Seal이 포함된 PBT Cap (내열 온도 : -40 ~ 200°C)
- 멸균시 반드시 뚜껑을 풀어 내압이 빠져갈 수 있도록 해야함
- LukeGL®

