

Flat bottom Vessel

평저 용기

Cat. No.	Model	Size (Ø × h)	Cap	EA Price ₩
R12-132-130	CLS-1437-B01	32 × 94 mm	GL32	448,200

- Addition and Sampling vessel with hose barbed 2-port drip tube assemblies
- GL32 규격의 Open top cap 포함
- Chemglass

Flat Bottom Vessel Only

Cat. No.	Model	Size (Ø × h)	EA Price ₩
R12-132-133	CLS-1437-B02	32 × 94 mm	76,500



Nessler Cylinder

비색 실린더, LukeGL®

Cat. No.	Model	용량	눈금 (ml)	Tolerance	EA Price ₩
C27-503-501	LK20201	50 ml	25 & 50	±0.4 ml	5,000
C27-503-503	LK20202	100 ml	50 & 100	±0.8 ml	9,200

- Borosilicate Glass 3.3
- Pouring Spout 를 가지고 있음
- LukeGL®



Nessler Tube with Stopper


공전 비색관, LukeGL®

Tall Form, 장형

Cat. No.	Model	용량	Size	EA Price ₩
T16-96-056	NT-T025S	25 ml	Ø13 × h250 mm	11,700
T16-96-062	NT-T050S	50 ml	Ø20 × h290 mm	14,600
T16-96-068	NT-T100S	100 ml	Ø24 × h370 mm	19,000

Short Form, 단형

Cat. No.	Model	용량	Size	EA Price ₩
T16-96-111	NT-S050S	50 ml	Ø24 × h205 mm	14,600
T16-96-119	NT-S100S	100 ml	Ø30 × h260 mm	19,000



APHA Nessler Color Comparison Tube, Kimble®

APHA 네슬러 비색관

- 재질 : 33 Expansion, Low Extractable Borosilicate Glass
- 하단에 그림자가 지지 않아 빛의 왜곡이 발생하지 않음
- 규격이 일치하는 6개가 한 세트임 (45325B 제품은 12개가 한 세트)
- DWK Lifescience / Kimble®

without Cap ①

Cat. No.	Model	Size (od × h)	눈금	PK Price ₩
T16-137-311	45315A-50	20 × 300 mm	50 ml	6ea 292,300
T16-137-315	45315A-100	24 × 375 mm	100 ml	6ea 403,300

with Cap ②

Cat. No.	Model	Size (od × h)	눈금	Joint	PK Price ₩
T16-137-319	45325A-50	20 × 312 mm	50 ml	19/10	6ea 539,300
T16-137-323	45325B-50100	24 × 390 mm	50 & 100 ml	24/12	12ea 1,570,000

Spare Cap ③

Cat. No.	Model	규격	EA Price ₩
T16-137-728	15182-1910	19/10	58,300
T16-137-731	15182-2412	24/12	78,600

