

5 mm Highest Quality NMR Tube, Kimble® 5 mm 정밀형 NMR 튜브

Cat. No.	Model	MHz	od (mm)	Concentricity	Camber	PK Price ₩
T16-137-351	897250-3000	900	4.97 mm (+0.0000 /-0.013), t 0.38 mm	5 µm	3.8 µm	5ea 400,000
T16-137-355	897245-3000	800				5ea 279,300
T16-137-359	897241-0000	600~700		13 µm	6 µm	5ea 218,900
T16-137-363	897240-0000	500				5ea 154,700
T16-137-367	897235-0000	400		25 µm	13 µm	5ea 113,300
T16-137-371	897230-0000	300		38 µm	25 µm	5ea 98,100
T16-137-375	897225-0000	200		51 µm	13 µm	5ea 98,100
T16-137-379	897220-0000	150		51 µm	25 µm	5ea 75,400
T16-137-387	897200-0000	100		76 µm	51 µm	5ea 52,900


- 재질 : Kimax® Borosilicate 3.3, ASTM Type I, Class A Glass
- 지름 : Ø5 mm, Length : 7" (178 mm)
- Ideal for Sealing Directly to Vacuum Manifold
- DWK Lifescience / Kimble®



Thrift Grade 5 mm NMR Tube, Kimble® 5 mm 경제형 NMR 튜브

Cat. No.	Model	MHz	Wall Thickness	Concentricity	Camber	PK Price ₩
T16-137-391	897150-0007	500	0.38 mm	13 µm	6 µm	5ea 71,700
T16-137-395	897140-0007	400				5ea 51,000
T16-137-399	897130-0007	300		38 µm	25 µm	5ea 44,500
T16-137-403	897120-0007	200		51 µm	13 µm	5ea 34,000
T16-137-407	897110-0007	>200			25 µm	5ea 31,000

- 재질 : Kimax® Borosilicate 5.1, ASTM Type I, Class B Glass
- 지름 : Ø5 mm, Length : 7" (178 mm)
- DWK Lifescience / Kimble®



5 mm Disposable Grade NMR Tube, Kimble® 5 mm 일회용 NMR 튜브

Cat. No.	Model	MHz	Wall Thickness	PK Price ₩
T16-137-415	897193-0000	100	0.38 mm	25ea 162,300
T16-137-419	897193-0050			50ea 184,900
T16-137-423	897193-7100			100ea 339,600

- 재질 : Kimax® Borosilicate 5.1, ASTM Type I, Class B Glass
- 지름 : Ø5 mm, Length : 7" (178 mm)
- PE Cap이 분리 포장되어져 있음
- DWK Lifescience / Kimble®



3 mm Highest Quality NMR Tube, Kimble® 3 mm 정밀형 NMR 튜브

Cat. No.	Model	MHz	Wall Thickness	Concentricity	Camber	PK Price ₩
T16-137-327	897840-0000	600~800	0.284 mm	0.0005"	0.00025"	5ea 192,400
T16-137-331	897835-0000	500		0.001"	0.0005"	5ea 166,000
T16-137-335	897830-0000	400		0.0015"	0.001"	5ea 192,700
T16-137-343	897805-0000	100		0.0025"	0.002"	5ea 158,200
T16-137-347	897800-0000	80		0.003"		5ea 174,800

- 재질 : Kimax® Borosilicate 3.3, ASTM Type I, Class B Glass
- 외경 : 3.0124 (+0.00 / -0.013) mm, 내경 : 2.4003 (+0.013 / +0.00) mm
- PE Cap 포함
- DWK Lifescience / Kimble®

