

5 mm NMR Tube for Routine NMR 5 mm NMR 튜브

- 재질 : N51A, ASTM Type I Class B Glass
- Samples are run under Room Temp with No Thermal Gradient
- Reaction Monitoring Scouting for a New Compound 및 Routine NMR용
- Auto-sampling 및 High-throughput NMR 튜브
- Length : 7" (178 mm)
- Chemglass / Norell

Precision NMR Tube

Cat. No.	Model	MHz	od (mm)	Concentricity	Camber	PK Price ₩	CS Price ₩
T16-132-336	C-505-P-7	200	4.97±0.030	±0.010 mm	±0.040 mm	5ea 51,600	5pk 226,600
T16-132-343	C-XR-55-7	300	4.97±0.025	±0.010 mm	±0.038 mm	5ea 53,300	5pk 253,200
T16-132-350	C-507-HP-7	400	4.97±0.013	±0.007 mm	±0.019 mm	5ea 68,200	-
T16-132-357	C-508-UP-7	500	4.97±0.013	±0.005 mm	±0.013 mm	5ea 95,200	-
T16-132-364	C-509-UP-7	600	4.97±0.006	±0.004 mm	±0.006 mm	5ea 198,800	-

High-throughput NMR Tube

Cat. No.	Model	MHz	od (mm)	Concentricity	Camber	PK Price ₩	CS Price ₩
T16-132-322	C-ST500-7	HTPLUS	4.97±0.070	±0.025mm	±0.075mm	5ea 21,400	20pk 407,500
T16-132-329	C-502-7	HT	4.97±0.050	±0.020mm	±0.070mm	5ea 23,600	10pk 224,200

- For Near Room Temp. Analysis of Routine Samples Exposed only to Slight Thermal Gradient

5 mm Amber NMR Tube 5 mm 갈색 NMR 튜브

Cat. No.	Model	MHz	od (mm)	Concentricity	Camber	EA Price ₩
T16-132-957	CS-5-500-7A	500	4.97±0.013	±0.0050 mm	±0.0130 mm	89,400

- 재질 : Borosilicate 3.3, ASTM Type I Class A Glass
- Typical Absorbance (ABS)
 - 0.5 ~ 1.0 at 650 nm ~ 300 nm
 - ≥ 1.0 at 300 nm ~ 190 nm
- 요청에 의해 다른 MHz 범위의 제품도 공급 가능함
- Chemglass / Norell

5 mm NMR Tube for High-Resolution NMR 고해상도 5 mm NMR 튜브

Cat. No.	Model	MHz	od (mm)	Concentricity	Camber	PK Price ₩
T16-132-481	CS-5-200-7	200	4.97±0.0025	±0.0090 mm	±0.0350 mm	5ea 64,100
T16-132-483	CS-5-300-7	300	4.97±0.0040	±0.0070 mm	±0.0250 mm	5ea 102,200
T16-132-485	CS-5-400-7	400	4.97±0.0050	±0.0070 mm	±0.0190 mm	5ea 132,700
T16-132-487	CS-5-500-7	500	4.97±0.0070	±0.0050 mm	±0.0130 mm	5ea 189,400
T16-132-489	CS-5-600-7	600	4.97±0.0070	±0.0040 mm	±0.006 mm	5ea 286,600
T16-132-491	CS-5-800-7	800	4.97±0.0090	±0.0025 mm	±0.0038 mm	5ea 416,700

- 재질 : Borosilicate 3.3, ASTM Type I Class A Glass
- 열충격에 강하고 팽창계수는 낮으므로 다양한 온도에서 사용
- Length : 7" (178 mm)
- Chemglass / Norell

5 mm Quartz NMR Tube 5 mm 석영 NMR 튜브

Cat. No.	Model	MHz	길이	PK Price ₩
T16-132-677	C-S-5-500-QTZ-7	500	178 mm (7")	5ea 646,300
T16-132-684	C-S-5-600-QTZ-7	600		5ea 817,500

- 재질 : Thin Wall Quartz, 99.8% SiO₂
- Boron NMR (0.1 ppm Boron) 용
- UV Catalyzed Reaction in the Region above 210 nm용
- Length : 7" (178 mm)
- Chemglass / Norell

