

B.O.D Bottle, 300 ml with Numbered & Un-numbered
300ml BOD 바틀, Robotic Stopper 포함

Cat. No.	Model	용량	Number	EA Price ₩	PK Price ₩
B08-503-941	BOD300RA	300 ml	Unnumbered	15,900	-
B08-503-943	BOD300RB		01 ~ 24	20,000	24ea 457,500
B08-503-945	BOD300RC		25 ~ 48	20,000	24ea 457,500
B08-503-947	BOD300RD		49 ~ 72	20,000	24ea 457,500
B08-503-949	BOD300RE		73 ~ 96	20,000	24ea 457,500
B08-503-951	BOD300RF		97 ~ 120	20,000	24ea 457,500
B08-503-953	BOD300RG		121 ~ 144	20,000	24ea 457,500
B08-503-955	BOD300RH		145 ~ 168	20,000	24ea 457,500
B08-503-957	BOD300RI		169 ~ 192	20,000	24ea 457,500
B08-503-959	BOD300RJ		193 ~ 216	20,000	24ea 457,500
B08-503-961	BOD300RK		217 ~ 240	20,000	24ea 457,500

- USP Type I, ASTM E438, Type I, Class A, Borosilicate Glass 3.3
- 5 Day Biochemical Oxygen Demand (BOD) Test 용
- Glass Robotic Stopper 기본 포함
- BOD 및 DO Meter의 Probe에 알맞은 Stopper Joint
- Size : Ø68×h140 mm, Neck id : Approx. 19 mm
- LukeGL®

Optional Accessory

Cat. No.	Model	용량	EA Price ₩
B08-503-901	ROS-BOD	Robotic Stopper, Spare	2,800
B08-503-903	CAP-BOD	Snap Fit Cap, PP for 300 ml BOD Bottle	700
B08-55-114	LK30012	Wire Rack for 300ml BOD Bottle, 12개용	40,200



B.O.D Bottle, 60 ml with Numbered & Un-numbered
60ml BOD 바틀, Robotic Stopper 포함

Cat. No.	Model	용량	Number	EA Price ₩	PK Price ₩
B08-503-911	BOD060RA	60 ml	Unnumbered	11,900	-
B08-503-913	BOD060RB		01 ~ 36	14,900	36ea 513,800
B08-503-915	BOD060RC		37 ~ 72	14,900	36ea 513,800
B08-503-917	BOD060RD		73 ~ 108	14,900	36ea 513,800

- USP Type I, ASTM E438, Type I, Class A, Borosilicate Glass 3.3
- 5 Day Biochemical Oxygen Demand (BOD) Test 용
- Glass Pennyhead Stopper 기본 포함
- BOD 및 DO Meter의 Probe에 알맞은 Stopper Joint
- Size : Ø43×h110 mm, Neck id : Approx. 26 mm
- LukeGL®

Optional Accessory

Cat. No.	Model	Description	EA Price ₩
B08-53-111	LK06020	Wire Rack for 60ml BOD Bottle, 20개용	55,700



Heavy Duty Glass Centrifuge Bottle, Kimble®
헤비월 유리 원심분리기용 바틀

Cat. No.	Model	용량	Size	EA Price ₩
B08-137-152	14720-200	200 ml	od60×h145 mm	105,300

- KIMAX® Borosilicate Glass 3.3, USP Type I, ASTM E438 Type I, Class A
- Small Bottom Area : Better Concentration, Max.RCF : 3550
- Pulp / Vinyl Coated Liner를 가진 Black Closure, GPI : 38-400
- DWK Lifescience / Kimble®

