


5 mm Valved NMR Tube for Vacuum & Reduced Pressure 5 mm 진공밸브형 NMR 튜브

Cat. No.	Model	MHz	id	Vacuum	EA Price ₩
T16-132-985	C-S5500-JY7	500	4.20 mm	138 ul/cm	482,000
T16-132-987	C-S5600-JY7	600			520,200


- 재질 : Borosilicate 3.3, ASTM Type I Class A Glass
- PTFE Valve & Fluoropolymer covered O-ring
- NMR Sample without Flame-sealing
- 사용가능조건 : (Vacuum) 10^{-7} kPa (10^{-6} torr) / (Pressure) 500 kPa (5 bar, 72 psi)
- Length : 7" (178 mm)
- Chemglass / Norell



5 mm NMR Tube with Screw Cap 5 mm 스크류캡 NMR 튜브

Cat. No.	Model	MHz	길이	EA Price ₩
T16-132-974	CS-5-500-SC-7	500	7" (178 mm)	139,300
T16-132-975	CS-5-600-SC-7	600		145,000


- 재질 : Borosilicate 3.3, ASTM Type I Class A Glass
- 8-425 Opentop Screwcap with PTFE/Silicon Septum
- 22 Gauge Standard Point Needle과 사용하기 적합
- Chemglass / Norell



PTFE Liner 테프론 라이너

Cat. No.	Model	Tube Size	Volume at 50 mm	EA Price ₩
T16-132-690	C-TL-5-7	5 mm	Approx. 0.35 ml	49,600
T16-132-697	C-TL-10-7	10 mm	Approx. 2.00 ml	89,300

- PTFE Plug Closure 와 같이 공급
- Hydrofluoric Acid, Ammonium Bifluoride 와 같은 Chemical Compound 와 같이 사용
- Chemglass / Norell



NMR Tube Cap NMR 튜브 캡

- 재질 : PE (Polyethylene)
- Chemglass / Norell

5 mm NMR Tube Cap

Cat. No.	Model	색상	PK Price ₩
T16-132-738	1842-01	Green	100ea 16,300
T16-132-742	1842-02	Red	100ea 16,300
T16-132-746	1842-03	Blue	100ea 16,300
T16-132-750	1842-04	Yellow	100ea 16,300
T16-132-754	1842-05	White	100ea 16,300
T16-132-757	1842-06	Black	100ea 16,300

5 mm NMR Tube Cap Package (in Box)

Cat. No.	Model	규격	PK Price ₩
T16-132-764	1842-10	6종류 컬러 5 mm용 캡	600ea 114,900

- Green, Red, Blue, Yellow, White, Black 각 100개씩 PP박스에 포장되어져 있음



PTFE NMR Tube Cap 테프론 NMR 튜브 캡

Cat. No.	Model	Tube Size	색상	PK Price ₩
T16-132-855	TC-5-PTFE	5 mm	White	50ea 553,300
T16-132-862	TC-10-PTFE	10 mm		50ea 835,500

- 재질 : PTFE
- Chemglass / Norell

