


Quartz Tube 석영 관

배송비 별도

| Cat. No. | Model | 외경 (mm) | 내경 (mm) | 두께 (mm) | 길이 (mm) | EA Price ₩ |
|------------|---------|-----------|-----------|----------|---------|------------|
| Q01-60-006 | QTS1001 | 3.2±0.2 | 1.7±0.2 | 0.75±0.2 | 1200 mm | 12,400 |
| Q01-60-012 | QTS1002 | 4.0±0.2 | 2.0±0.2 | 1.0±0.2 | | 12,400 |
| Q01-60-018 | QTS1003 | 5.0±0.2 | 3.0±0.2 | 1.0±0.2 | | 12,400 |
| Q01-60-024 | QTS1004 | 6.0±0.2 | 4.0±0.2 | 1.0±0.2 | | 12,400 |
| Q01-60-030 | QTS1005 | 7.0±0.2 | 5.0±0.2 | 1.0±0.2 | | 12,400 |
| Q01-60-036 | QTS1006 | 8.0±0.2 | 5.0±0.2 | 1.5±0.2 | | 12,900 |
| Q01-60-039 | QTS1106 | 8.0±0.2 | 6.0±0.2 | 1.0±0.2 | | 12,400 |
| Q01-60-042 | QTS1007 | 9.0±0.2 | 7.0±0.2 | 1.0±0.2 | | 12,400 |
| Q01-60-048 | QTS1008 | 10.0±0.2 | 8.0±0.2 | 1.0±0.2 | | 12,400 |
| Q01-60-054 | QTS1009 | 11.0±0.3 | 9.0±0.3 | 1.0±0.3 | | 12,400 |
| Q01-60-067 | QTS1109 | 12.0±0.3 | 8.0±0.3 | 2.0±0.3 | | 18,100 |
| Q01-60-060 | QTS1010 | 12.0±0.3 | 10.0±0.3 | 1.0±0.3 | | 13,100 |
| Q01-60-066 | QTS1011 | 13.0±0.3 | 11.0±0.3 | 1.0±0.3 | | 14,300 |
| Q01-60-072 | QTS1012 | 14.0±0.3 | 12.0±0.3 | 1.0±0.3 | | 14,300 |
| Q01-60-074 | QTS1112 | 15.0±0.3 | 11.0±0.3 | 1.5±0.3 | | 21,300 |
| Q01-60-075 | QTS1113 | 15.0±0.3 | 12.0±0.3 | 1.5±0.3 | | 18,200 |
| Q01-60-078 | QTS1013 | 15.0±0.3 | 13.0±0.3 | 1.0±0.3 | | 14,900 |
| Q01-60-081 | QTS1114 | 16.0±0.3 | 12.0±0.3 | 2.0±0.3 | | 22,500 |
| Q01-60-084 | QTS1014 | 16.0±0.3 | 14.0±0.3 | 1.0±0.3 | | 16,100 |
| Q01-60-090 | QTS1015 | 17.0±0.3 | 15.0±0.3 | 1.0±0.3 | | 17,200 |
| Q01-60-093 | QTS1115 | 18.0±0.3 | 15.0±0.3 | 1.5±0.3 | | 21,300 |
| Q01-60-096 | QTS1016 | 18.0±0.3 | 16.0±0.3 | 1.0±0.3 | | 18,200 |
| Q01-60-102 | QTS1017 | 19.0±0.3 | 17.0±0.3 | 1.0±0.3 | | 19,300 |
| Q01-60-105 | QTS1118 | 20.0±0.3 | 16.0±0.3 | 2.0±0.3 | | 27,600 |
| Q01-60-108 | QTS1018 | 20.0±0.3 | 17.0±0.3 | 1.5±0.3 | | 22,800 |
| Q01-60-114 | QTS1019 | 20.5±0.3 | 18.0±0.3 | 1.25±0.3 | | 22,000 |
| Q01-60-120 | QTS1020 | 22.0±0.3 | 19.0±0.3 | 1.5±0.3 | | 25,100 |
| Q01-60-126 | QTS1021 | 23.0±0.3 | 20.0±0.3 | 1.5±0.3 | | 26,300 |
| Q01-60-132 | QTS1022 | 24.0±0.3 | 21.0±0.3 | 1.5±0.3 | | 27,600 |
| Q01-60-135 | QTS1122 | 25.0±0.3 | 21.0±0.3 | 2.0±0.3 | | 31,400 |
| Q01-60-138 | QTS1023 | 25.0±0.3 | 22.0±0.3 | 1.5±0.3 | | 29,000 |
| Q01-60-144 | QTS1024 | 26.0±0.3 | 23.0±0.3 | 1.5±0.3 | | 34,500 |
| Q01-60-156 | QTS1026 | 28.0±0.3 | 25.0±0.3 | 1.5±0.3 | 37,600 | |
| Q01-60-168 | QTS1028 | 30.0±0.3 | 27.0±0.3 | 1.5±0.3 | 35,800 | |
| Q01-60-174 | QTS1029 | 32.0±0.5 | 28.0±0.5 | 2.0±0.5 | 44,500 | |
| Q01-60-180 | QTS1030 | 33.0±0.5 | 29.0±0.5 | 2.0±0.5 | 52,000 | |
| Q01-60-186 | QTS1031 | 34.0±0.5 | 30.0±0.5 | 2.0±0.5 | 63,400 | |
| Q01-60-192 | QTS1032 | 35.0±0.5 | 31.0±0.5 | 2.0±0.5 | 63,400 | |
| Q01-60-198 | QTS1033 | 36.0±0.5 | 32.0±0.5 | 2.0±0.5 | 69,600 | |
| Q01-60-210 | QTS1035 | 38.0±0.5 | 34.0±0.5 | 2.0±0.5 | 79,600 | |
| Q01-60-228 | QTS1038 | 40.0±0.5 | 36.0±0.5 | 2.0±0.5 | 85,900 | |
| Q01-60-240 | QTS1040 | 42.0±0.75 | 38.0±0.75 | 2.0±0.75 | 101,100 | |
| Q01-60-252 | QTS1042 | 44.0±0.75 | 40.0±0.75 | 2.0±0.75 | 117,300 | |
| Q01-60-258 | QTS1043 | 45.0±0.75 | 41.0±0.75 | 2.0±0.75 | 110,500 | |
| Q01-60-264 | QTS1044 | 48.0±0.75 | 44.0±0.75 | 2.0±0.75 | 182,000 | |
| Q01-60-270 | QTS1045 | 50.0±1.0 | 46.0±1.0 | 2.0±1.0 | 169,000 | |
| Q01-60-276 | QTS1046 | 52.0±1.0 | 48.0±1.0 | 2.0±1.0 | 195,000 | |
| Q01-60-282 | QTS1047 | 55.0±1.0 | 50.0±1.0 | 2.5±1.0 | 195,000 | |
| Q01-60-288 | QTS1048 | 60.0±1.0 | 55.0±1.0 | 2.5±1.0 | 227,500 | |
| Q01-60-294 | QTS1049 | 65.0±1.0 | 60.0±1.0 | 2.5±1.0 | 253,500 | |
| Q01-60-300 | QTS1050 | 70.0±1.5 | 65.0±1.5 | 2.5±1.5 | 289,200 | |
| Q01-60-303 | QTS1150 | 75.0±1.5 | 70.0±1.5 | 2.5±1.5 | 325,000 | |
| Q01-60-306 | QTS1051 | 80.0±1.5 | 74.0±1.5 | 3.0±1.5 | 377,000 | |
| Q01-60-309 | QTS1151 | 90.0±1.5 | 84.0±1.5 | 3.0±1.5 | 445,400 | |
| Q01-60-312 | QTS1052 | 100.0±2.0 | 94.0±2.0 | 3.0±0.3 | 611,000 | |

- 1250℃ 내열, 상용 온도 1,100℃
- 높은 기계적 강도와 화학적 내구성, 낮은 열 팽창율
- 강산과 약염기에서 매우 안정 (강염기에서는 불안정)
- 넓은 온도 범위에서도 정확한 양을 측정할 수 있기 때문에 특히 실험기구에 많이 사용
- LK Labkorea



T

Tube
Quartz