

Tubing for Peristaltic Pump, Masterflex®
정량펌프용 튜빙, Platinum Silicon

Cat. No.	Model	id (mm)	Hose Barb (mm)	Max. Pres.(Mpa)		EA Price ₩
				연속	불연속	
T23-55-011	96410-13	0.8	1.6	0.17	0.27	122,400
T23-55-015	96410-14	1.6	1.6			150,000
T23-55-019	96410-16	3.1	3.2			144,600
T23-55-023	96410-25	4.8	4.8	0.14	0.24	189,500
T23-55-031	96410-17	6.4	6.4	0.1	0.14	200,500
T23-55-039	96410-18	7.9	9.5	0.07	0.1	234,500
T23-55-527	96410-15	4.8	4.8	0.17	0.27	196,400
T23-55-535	96410-24	6.4	6.4			267,200
T23-55-543	96410-35	7.9	9.5	0.14	0.24	317,600
T23-55-551	96410-36	9.7	9.5	0.1	0.14	362,600

- 사용 가능 온도 : -60 ~ 200°C
- 길이 : 7.6 m
- 멸균 가능, Smooth inner surface, minimizes particle entrapment
- Do not use with concentrated acid and bases, organic solvents and oils
- 정량펌프와 같이 사용하기 위한 Tubing
- 보다 다양하고 자세한 정량펌프용 Tubing은 Pump란 참조
- **Masterflex®**



- Masterflex® 정량 펌프용 L/S Tubing은 820 Page 참조
- Masterflex® 정량 펌프용 I/P Tubing은 823 Page 참조

Viton Tubing
바이톤 튜빙

Cat. No.	Model	id (mm)	od (mm)	Max. psi at 70°C	EA Price ₩
T23-55-911	643401	1.6	3.2	23	131,300
T23-55-915	643402	3.2	4.8	21	121,800
T23-55-919	643403	4.8	6.4	19	183,800
T23-55-923	643404	6.4	8.0	17	218,800
T23-55-927	643405	8.0	11.2	20	456,000
T23-55-931	643406	9.6	12.8	18	538,300
T23-55-935	643407	12.8	16.0	16	524,400
T23-55-939	643408	16.0	22.4	17	1,341,700
T23-55-943	643409	19.1	25.4	15	1,494,000

- Viton® 재질의 검정색 튜빙
- Acid나 Solvent와 같이 사용
- 사용 온도 범위 : -32 ~ 204°C
- 7.6 m/box
- LK Labkorea



Triple Insulated Flexible Tubing
고온용 튜빙, Stainless Steel

Cat. No.	Model	Female End	길이	EA Price ₩
T23-132-101	CG-1969-B-01	M16	30 cm	715,200
T23-132-103	CG-1969-B-02		60 cm	748,600
T23-132-105	CG-1969-B-04		1.2 m	1,064,500
T23-132-107	CG-1969-B-06		1.8 m	1,314,000
T23-132-109	CG-1969-B-08		2.4 m	1,746,900

- 사용온도 : -80 ~ 200°C
- Triple Insulated Flexible Stainless Steel Tubing
- **Chemglass**

