


PTFE Tubing

PTFE 테프론 튜빙, 불투명

- 백색의 불투명하고 단단한 내화학, 내열성 튜브
- -200°C + 280°C 내열
- LK Labkorea

Metric PTFE 튜빙

Cat. No.	Model	id	od	길이	EA Price ₩
T23-140-304	TF-0101	Ø1 mm	Ø2 mm	10 m	29,000
T23-140-308	TF-0102	Ø1.5 mm	Ø2.5 mm		25,500
T23-140-312	TF-0203	Ø2 mm	Ø3 mm		25,500
T23-140-316	TF-0304	Ø3 mm	Ø4 mm		25,500
T23-140-320	TF-0406	Ø4 mm	Ø6 mm		34,700
T23-140-324	TF-0507	Ø5 mm	Ø7 mm		46,500
T23-140-328	TF-0608	Ø6 mm	Ø8 mm		50,200
T23-140-332	TF-0709	Ø7 mm	Ø9 mm		65,100
T23-140-336	TF-0810	Ø8 mm	Ø10 mm		60,300
T23-140-340	TF-0911	Ø9 mm	Ø11 mm		66,100
T23-140-344	TF-1012	Ø10 mm	Ø12 mm		71,900
T23-140-348	TF-1214	Ø12 mm	Ø14 mm		82,700
T23-140-350	TF-1416	Ø14 mm	Ø16 mm		112,100
T23-140-352	TF-1517	Ø15 mm	Ø17 mm		106,100
T23-140-356	TF-1618	Ø16 mm	Ø18 mm		115,700
T23-140-360	TF-1820	Ø18 mm	Ø20 mm		127,800



Inch-type PTFE 튜빙

Cat. No.	Model	id	od	길이	EA Price ₩
T23-140-370	TF-IN032	1/32"	1/16"	10 m	46,300
T23-140-374	TF-IN070	1/16"	1/8"		76,400
T23-140-378	TF-IN155	5/32"	1/4"		90,200
T23-140-380	TF-IN175	4.35 mm	1/4"		81,000
T23-140-382	TF-IN195	1/4"	3/8"		165,400
T23-140-384	TF-IN199	7.53 mm	3/8"		140,000
T23-140-386	TF-IN215	3/8"	1/2"		216,300
T23-140-390	TF-IN245	17mm	3/4"		200,400
T23-140-394	TF-IN295	22.6mm	1"		336,800


PFA Tubing

PFA 테프론 튜빙, 투명

- 투명하고 단단한 내화학, 내열성 튜브
- -200°C + 260°C 내열
- LK Labkorea

Metric PFA 튜빙

Cat. No.	Model	id	od	길이	EA Price ₩
T23-140-416	PF-0203	2 mm	3 mm	10 m	72,300
T23-140-420	PF-0204	2 mm	4 mm		123,400
T23-140-424	PF-0304	3 mm	4 mm		81,300
T23-140-428	PF-0406	4 mm	6 mm		136,500
T23-140-432	PF-0608	6 mm	8 mm		226,700
T23-140-436	PF-0810	8 mm	10 mm		284,600
T23-140-440	PF-1012	10 mm	12 mm		444,200
T23-140-444	PF-1214	12 mm	14 mm		511,300
T23-140-448	PF-1416	14 mm	16 mm		583,800



Inch-type PFA 튜빙

Cat. No.	Model	id	od	길이	EA Price ₩
T23-140-470	PF-IN070	1/16"	1/8"	10 m	115,700
T23-140-472	PF-IN162	1/8"	1/4"		208,800
T23-140-474	PF-IN160	5/32"	1/4"		196,700
T23-140-476	PF-IN165	4.35 mm	1/4"		173,500
T23-140-478	PF-IN195	1/4"	3/8"		393,300
T23-140-480	PF-IN199	7.53 mm	3/8"		283,200
T23-140-482	PF-IN215	3/8"	1/2"		548,300
T23-140-486	PF-IN245	17 mm	3/4"		587,400
T23-140-488	PF-IN251	23 mm	1"		705,700