

Screw Cap & Pouring Ring 랩 바틀용 스크류 캡 & 링

- PP 재질의 기본형과 PBT 재질의 고온용 두가지임
- 내열온도 : -40 ~ 140°C (PP), -100°C ~ 200°C (PBT)
- KAVALIERGLASS

PP Screw Cap & Pouring Ring, Blue

Cat. No.	Model	Type	Size	EA Price ₩
B08-985-500	9180002279	Screw Cap Blue ①	GL32	3,200
B08-985-501	9180001452		GL45	1,300
B08-985-503	9180002271		GL80	2,800
B08-985-509	9180002280	Pouring Ring Blue ②	GL32	700
B08-985-511	9180001449		GL45	600
B08-985-513	9180002272		GL80	1,300



PBT Screw Cap & Pouring Ring, Red

Cat. No.	Model	Description	EA Price ₩
B08-985-601	9180001451	GL45 Screw Cap, Red ③	2,700
B08-985-611	9180001447	GL45 Pouring Ring, Red ④	2,500

- 9180001451 : PTFE / Silicon Gasket과 같이 사용 가능함 (불포함)

PTFE/Silicon Gasket for GL45 ⑤

Cat. No.	Model	Description	EA Price ₩
J02-460-139	H973-45	PTFE/Silicon Gasket, GL45용	4,700



GL45 Bottle Cap, Kimble® GL45 바틀 캡

- Autoclavable
- DWK Lifescience / Kimble®

PP Screw Thread Cap ①

Cat. No.	Model	색상	EA Price ₩	PK Price ₩
B08-137-530	14395C-45	Blue	6,200	10ea 58,000
B08-137-533	14395C-451	White	7,000	10ea 65,100
B08-137-536	14395C-453	Orange	7,000	10ea 65,100
B08-137-539	14395C-455	Gray	7,000	10ea 66,100
B08-137-542	14395C-457	Green	7,000	10ea 65,100
B08-137-545	14395C-459	Pink	7,000	10ea 65,100

- Internal Molded Seal Ring

PP Pour Ring ②

Cat. No.	Model	색상	EA Price ₩	PK Price ₩
B08-137-527	14395P-45	Clear	2,600	10ea 24,300

- Max. Temp. : 140°C



Adapter for HPLC Reservoir Bottle HPLC 리저버 바틀용 어댑터

Cat. No.	Model	재질	Description	EA Price ₩
B08-137-224	953905-0000	PTFE	GL45 to 40/35 ①	506,400
B08-137-228	953907-0000	PBT	GL45 to Solvent Bottle ②	279,500
B08-137-554	953908-0000	PTFE	Sealing Cap for GL45 ③	379,200

- Autoclavable PTFE & PBT
- 953905-0000 : GL45 Bottle를 Vacuum Filter Holder의 40/35 Male Joint와 연결 가능
- 953907-0000 : Solvent Bottle Adapter
- 953908-0000 : PTFE GL45 HPCL Storage Cap with O-ring for HPCL Reservoirs
- DWK Lifescience / Kimble®

