

UV Handy Lamp, LK Lab®

자외선 핸드 램프

- Single 또는 Dual 파장의 UV Lamp
- Intensity는 15 cm 거리에서의 측정 값
- 모든 램프는 필터링 되며 형광 검출에 사용하기 알맞음
- Unfiltered Lamp, 동일파장(x2) 주문제작 가능
- LK Labkorea

4W UV Lamp



Cat. No.	Model	Tubes (Watt)	Wavelength (nm)	Intensity (uW/cm ²)	EA Price ₩
U01-02-101	LT-4L	1×4 W	365	400	397,600
U01-02-105	LT-4M		312	400	409,600
U01-02-109	LT-4C		254	340	397,600
U01-02-121	LT-4LC	2×4 W	365/254	400/340	421,700
U01-02-125	LT-4LM		365/312	400/400	433,700
U01-02-129	LT-4MC		312/254	400/340	433,700

6W UV Lamp

Cat. No.	Model	Tubes (Watt)	Wavelength (nm)	Intensity (uW/cm ²)	EA Price ₩
U01-02-141	LT-6L	1×6 W	365	610	421,700
U01-02-145	LT-6M		312	580	433,700
U01-02-149	LT-6C		254	400	421,700
U01-02-161	LT-6LC	2×6 W	365/254	610/400	445,800
U01-02-165	LT-6LM		365/312	610/580	469,900
U01-02-169	LT-6MC		312/254	580/400	469,900

15W UV Lamp

Cat. No.	Model	Tubes (Watt)	Wavelength (nm)	Intensity (uW/cm ²)	EA Price ₩
U01-02-181	LT-115L	1×15 W	365	1350	783,100
U01-02-185	LT-115M		312	1100	843,400
U01-02-189	LT-115C		254	930	783,100
U01-02-201	LT-215LC	2×15 W	365/254	1350/1100	903,600
U01-02-205	LT-215LM		365/312	1350/1800	933,700
U01-02-209	LT-215MC		312/254	1800/930	933,700
U01-02-221	LT-215L	2×15 W	365	2300	903,600
U01-02-225	LT-215M		312	3000	963,900
U01-02-229	LT-215C		254	1780	903,600



UV Ultraviolet Lamp

자외선 램프

- 같은 용량의 UV Lamp의 Spare Lamp
- 조명용 형광등 기구와 같이 사용 가능함
- LK Labkorea

UV-C (254 nm), 살균 램프

Cat. No.	Model	Lamp	Size (Ø×L, mm)	UV Output	EA Price ₩
U01-133-264	T-4C	4 W	15.5×134.5	0.8	21,700
U01-133-285	T-6C	6 W	15.5×210.5	1.7	31,400
U01-133-290	T-8C	8 W	15.5×287	2.5	62,700
U01-133-296	T-10C	10 W	25.5×330	2.7	14,500
U01-133-306	T-15C	15 W	25.5×436	4.9	63,300
U01-133-322	T-20C	20 W	32.5×580	7.5	94,000
U01-133-340	T-30C	30 W	25.5×893	13.4	101,900
U01-133-365	T-40C	40 W	32.5×1198	19.8	109,700

- 자외선 254 nm의 파장역 (UV-C)을 방사
- DNA를 파괴하는 자외선, 세균을 사멸 시키는 작용을 함

UV-B (312 nm)

Cat. No.	Model	Lamp	Size (Ø×L, mm)	UV Output	EA Price ₩
U01-133-464	T-4M	4 W	15.5×134.5	0.8	108,500
U01-133-470	T-6M	6 W	15.5×210.5	1.7	77,100
U01-133-478	T-15M	15 W	25.5×436	4.9	94,800



UV-A (365 nm)

Cat. No.	Model	Lamp	Size (Ø×L, mm)	UV Output	EA Price ₩
U01-133-488	T-4L	4 W	15.5×134.5	0.8	31,400
U01-133-496	T-6L	6 W	15.5×210.5	1.7	39,200
U01-133-500	T-8L	8 W	15.5×287	2.5	62,700
U01-133-504	T-15L	15 W	25.5×436	4.9	77,100