

## Handheld Multi-Range Vacuum Gauge 휴대용 진공 게이지

Cat. No.	Model	Vacuum Range	Connection	EA Price ₩
V03-750-932	PIZA111	1050 ~ 10 <sup>-3</sup> mbar	KF16 Flange	2,067,500

- Combination Sensor : Piezo and Pirani
- Resolution : 1 mbar (1050 ~ 100), 0.5 mbar (99 ~ 1), 0.1 mbar (0.9 ~ 0.1), 0.01 mbar (0.09 ~ 0.01), 0.001 mbar (0.009 ~ 0.001)
- Accuracy : < 0.8% FS, Reproducibility : < 0.5%
- Max. Pressure : 2bar, Size : w80×d36×h125 mm
- Welch-Ilmvac



## Thermocouple Vacuum Gauge 온도 보정식 진공 게이지, Analog Type

Cat. No.	Model	Description	EA Price ₩
V08-132-231	AF-0301-01	Analog Gauge, Panel Mount	1,394,800
V08-132-235	AF-0301-02	Analog Gauge, Benchtop	2,165,900

- 대기상태에서 10<sup>-3</sup> Torr 영역까지 표시되며, 가장 널리 사용되고 있는 모델
- AF-0302-02 Thermocouple(Sensor) 기본 포함, 8-Pin, Hastings Style
- AF-0306 Adapter와 같이 사용가능 (1209 page 참조)
- Chemglass

Specification	
Range	1 × 10 <sup>-3</sup> torr to 2 torr
Recorder Output	0 ~ 11 mV (non - linear, Monotonic)
Lowest Indication	1 × 10 <sup>-3</sup> torr
외부 사용 온도 조건	10 ~ 40°C
Transducer Cable	10 feet (3.04 m) attached
Power Cable	6 feet (1.83 m) supplied
Net Weight	14 oz. (0.40 kg)
Power	220 V
Power	90 ~ 264V, 50/60 Hz



AF-0301-01



AF-0301-02

### Spare Accessory

Cat. No.	Model	Description	EA Price ₩
V08-132-500	AF-0302-02	Thermocouple Only, 1/8" Male NPT	275,100

## Analog Manometer 조절식 마노메타

Cat. No.	Model	mbar Range	Control Valve	EA Price ₩
V03-53-122	3000	0 ~ 1020	O (진공 조절식)	525,300
V03-53-132	3001		X	352,200

- 수은을 전혀 사용하지 않은 마노메타
- Tube Connection : OD 8 mm
- Size : w80×d80×h150 mm
- Van Der Heijden



## Vacuum Gauge 압력/진공 게이지

Cat. No.	Model	Measuring Range	EA Price ₩
V03-83-222	5802000	2 ~ 180 mbar	228,700

- Absolute Pressure 측정용
- Hard-wood Stand 위에 Adjustable Mirror Scale
- Mercury(수은) 포함
- LK Labkorea

