

## Double Bank Vacuum Manifold, LukeGL® 이중 메니폴드, 테프론 니들 밸브형

- Borosilicate Glass 3.3
- J-Young High Vacuum Precision PTFE Needle Valve : 10<sup>-6</sup> Torr
- 진공라인 및 Gas 투입라인
- 요청에 의해 규격외 사이즈도 공급 가능
- LukeGL®

### Two End Open (Cross) ①


Cat. No.	Model	Port 수	Connection	EA Price ₩
V03-132-580	4435-21T	3	Hose Nipple	867,500
V03-132-586	4435-22T	4		1,180,700
V03-132-592	4435-23T	5		1,566,300

### Three End Open ②

Cat. No.	Model	Port 수	Connection	EA Price ₩
V03-132-640	4435-31T	3	Hose Nipple	867,500
V03-132-646	4435-32T	4		1,180,700
V03-132-652	4435-33T	5		1,566,300

### Four End Open ③

Cat. No.	Model	Port 수	Connection	EA Price ₩
V03-132-676	4435-41T	3	Hose Nipple	867,500
V03-132-682	4435-42T	4		1,180,700
V03-132-688	4435-43T	5		1,566,300



## Double Bank Vacuum Manifold with Glass Stopcock 이중 메니폴드, 유리콕크형, LukeGL®

- Borosilicate Glass 3.3
- 진공라인 및 Inert Gas 투입라인
- Glass Stopcock : Bore 4.0 mm
- 요청에 의해 규격외 Size 공급 가능
- LukeGL®

### Two End Open (Cross) - Hose Nipple ④

Cat. No.	Model	Port 수	Connection	EA Price ₩
V03-10-062	L444011	3	Hose Nipple	408,600
V03-10-068	L444012	4		528,200
V03-10-074	L444013	5		663,700

### Three End Open - Hose Nipple ⑤

Cat. No.	Model	Port 수	Connection	EA Price ₩
V03-10-122	L444101	3	Hose Nipple	408,600
V03-10-128	L444102	4		528,200
V03-10-134	L444103	5		663,700

### Four End Open - Hose Nipple ⑥

Cat. No.	Model	Port 수	Connection	EA Price ₩
V03-10-158	L444141	3	Hose Nipple	413,500
V03-10-164	L444142	4		528,200
V03-10-170	L444143	5		663,700

### Four End Open - Standard Joint 24/40 ⑦

Cat. No.	Model	Port 수	Connection	EA Price ₩
V03-10-181	L444161	3	Socket 24/40	481,900
V03-10-183	L444162	4		626,500
V03-10-185	L444163	5		777,100

