



**Vacuum Manifold, AIRFREEE®**
**진공 메니폴드, 5-port**



Cat. No.	Model	Description	EA Price ₩
V08-132-301	AF-0450-01	Right Hand	11,576,600
V08-132-305	AF-0450-02	Left Hand	11,576,600

- Vacuum Bank : ID 40 mm, 34/25 Ball Joint
- 전체 사이즈 : w1016×h610 mm
- 구성
  - A : 4 mm Solid Double Oblique High Vacuum Stopcock (Pressure Retaining Device, 24/40 Inner Joint)
  - B : 35/25 Socket Joint Cap (Vacuum Bank 마개)
  - C : Valve Gauging Port (12/5 Ball Joint)
  - D : Inert Gas Line (Gas Source 연결부)
  - E : Valve, Cross Over Line (Evacuation of the entire system)
  - F : 14/20 Outer Joint Cap (Gas Line 마개)
  - G : Pre-trap (od51×h250 mm, 상/하부 연결부 28/15 Ball Joint)
  - H : Main Trap Body (od58×h235 mm, 상/하부 연결부 75/50 Ball Joint)
  - I : Valved Gauging port (Pre-trap & Main Trap 연결부 12/5 Ball Joint)
  - J : Pump 연결부 (od1"Tube)      - K : Isolation Valve (0 ~ 12 mm)
- Complete Set
  - 5 Port Vacuum Manifold (1ea), Pre-trap (1ea), Main Trap (1ea), #28 Pinch Clamp (2ea), #35 & #75 Pinch Clamp (각1ea), 1/2" Spring (2ea), 14/20 Outer Stopper (1ea), 34/25 Socket Stopper (1ea)
- Chemglass


**Vacuum Manifold, AIRFREEE®**
**진공 메니폴드, 4-port**

Cat. No.	Model	Description	EA Price ₩
V08-132-311	AF-0451-01	Vacuum Manifold, Right Hand	10,007,900
V08-132-315	AF-0451-02	Vacuum Manifold, Left Hand	10,007,900

- 5-port Line과 유사함      • 전체 사이즈 : w864×h610 mm
- 구성 차이점
  - A : 0 ~ 4 mm Valve      - B : #30 O-ring Joint Cap (Vacuum Bank 마개)
  - C : #9 O-ring Gauging Port      - D : #9 O-ring Joint Cap (Gas Line 마개)
- Complete Set
  - 4 Port Vacuum Manifold (1ea), Pre-trap (1ea), Main Trap (1ea), #28 Pinch Clamp (2ea), #35 Pinch Clamp (1ea), #50 Pinch Clamp (1ea), #75 Pinch Clamp (1ea), #18 Pinch Clamp (4ea), #112 Viton O-ring (4ea), #9 & #30 O-ring Joint Stopper (각1ea), #223 Viton O-ring (1ea), #9 O-ring Joint to 3/8" Tubulation Adapter (1ea)
- Chemglass


**Spare Part for Vacuum Manifold, AIRFREEE®**
**공통 부품 약세사리, Chemglass**

Cat. No.	Model	Description	EA Price ₩
V08-132-412	AF-0450-03	28/15 Pre-trap Only	341,600
V08-132-416	AF-0450-04	75/50 Main Trap Only	258,100
V08-132-422	AF-0451-05	#9 O-ring Cap Stopper	29,400
V08-132-426	AF-0451-06	#30 O-ring Cap Stopper	61,400
V08-132-429	AF-0451-08	#9 O-ring Joint to 9 mm Tubulation	38,900
V08-132-905	CG-305-112	#112 Viton O-ring, (12/pk)	21,400
V08-132-912	CG-305-223	#223 Viton O-ring, (3/pk)	18,600

Cat. No.	Model	Description	EA Price ₩
V08-132-620	CG-150-21	#75 SUS Pinch Clamp	136,300
V08-132-626	CG-150-08	#50 SUS Pinch Clamp	73,600
V08-132-623	CG-150-06	#35 SUS Pinch Clamp	30,400
V08-132-630	CG-150-05	#28 SUS Pinch Clamp	24,900
V08-132-634	CG-150-04	#18 SUS Pinch Clamp	18,400
V08-132-582	CG-961-01	0 ~ 4 mm Stem & Control Knob	74,300
V08-132-968	CG-961-03	0 ~ 12 mm Stem & Control Knob	84,000
V08-132-981	CG-110-04	1/2" SUS Spring (12/pk)	40,400

- Diffusion Pump 1210 Page 참조
- Dewar Flask 315 Page 참조

