

Sealing Ring for GL Screw Cap GL 스크류캡용 실링 링


- Silicon Ring with PTFE Washer
- One-side or Double-side와 Inner Hole에 PTFE 처리된 실링 링
- BOLA

One-sided PTFE Gasket


Cat. No.	Model	규격	Bore (mm)	EA Price ₩	PK Price ₩
J02-460-051	H975-02	GL14	5.5 ~ 6.5	2,300	10ea 22,900
J02-460-055	H975-04	GL18		6,400	10ea 60,800
J02-460-059	H975-06		7.5 ~ 9.0	5,800	10ea 54,900
J02-460-063	H975-10		9.0 ~ 11.0	5,100	10ea 47,700
J02-460-067	H975-12	GL25	7.5 ~ 9.0	7,100	10ea 68,200
J02-460-071	H975-14		9.0 ~ 11.0	6,600	10ea 63,300
J02-460-075	H975-18	GL32	11.0 ~ 13.0	7,700	10ea 74,100
J02-460-079	H975-20		9.0 ~ 11.0	11,700	10ea 111,800
J02-460-083	H975-22		11.0 ~ 13.0	10,200	10ea 97,300
J02-460-087	H975-26	GL45	13.0 ~ 15.0	9,800	10ea 93,300
J02-460-091	H975-28		15.0 ~ 17.0	9,800	10ea 93,300
J02-460-095	H975-30	GL45	17.0 ~ 19.0	11,800	10ea 112,000
J02-460-099	H975-34		25.0 ~ 27.0	19,300	10ea 183,700
J02-460-103	H975-36		31.0 ~ 33.0	16,600	10ea 157,500

Double-sided PTFE Gasket

Cat. No.	Model	규격	Bore (mm)	EA Price ₩	PK Price ₩
J02-460-011	H977-08	GL14	5.5 ~ 6.5	10,700	10ea 101,300
J02-460-015	H977-16	GL18	5.5 ~ 6.5	18,000	5ea 85,300
J02-460-019	H977-18		7.5 ~ 9.0	14,500	5ea 68,700
J02-460-023	H977-20		9.0 ~ 11.0	8,000	10ea 75,200
J02-460-027	H977-28	GL25	7.5 ~ 9.0	14,300	5ea 68,100
J02-460-031	H977-32		9.0 ~ 11.0	16,300	5ea 77,600
J02-460-035	H977-36	GL25	11.0 ~ 13.0	22,500	5ea 107,100



One Side




Double Side

PTFE/Silicon Gasket for GL Screw Cap GL 스크류 캡용 가스켓

- 재질 : PTFE, Silicon-elastomer
- 한쪽면이 PTFE 처리된 GL 스크류캡용 가스켓
- BOLA

Cat. No.	Model	규격	규격	EA Price ₩
J02-460-131	H973-14	GL 14	Ø13.0 x t3.3 mm	1,000
J02-460-133	H973-18	GL 18	Ø16.8 x t3.3 mm	1,700
J02-460-135	H973-25	GL 25	Ø23.5 x t3.3 mm	3,500
J02-460-137	H973-32	GL 32	Ø30.2 x t3.3 mm	3,900
J02-460-139	H973-45	GL 45	Ø43.2 x t3.3 mm	4,700



Silicon Rubber Septa 실리콘 고무 셉타

- 재질 : Silicon
- LK Labkorea

Cat. No.	Model	규격	두께	EA Price ₩
J02-83-062	1341214	GL14	2 mm	800
J02-83-065	1341618	GL18		800
J02-83-068	1342225	GL25		1,100
J02-83-071	1342932	GL32		1,400
J02-83-074	1344245	GL45	3 mm	1,500

