

## Jointed Vacuum Trap 진공 트랩, LukeGL®

Cat. No.	Model	Body Size	Joint	EA Price ₩
V03-132-261	VT-J024	∅28×h200 mm	24/40	38,200
V03-132-267	VT-J029	∅33×h220 mm	29/32	48,200
V03-132-273	VT-J034	∅40×h240 mm	34/45	49,400
V03-132-279	VT-J045	∅50×h280 mm	45/50	62,200

- Borosilicate Glass 3.3
- 분리형 진공 트랩으로 응축액 제거 용이
- LukeGL®



## Vacuum Trap 일체형 진공 트랩, LukeGL®

- Borosilicate Glass 3.3
- LukeGL®

### 기본형

Cat. No.	Model	Size	EA Price ₩
V03-132-202	4510-01	∅25×h200 mm	24,800
V03-132-208	4510-02	∅32×h250 mm	28,700
V03-132-214	4510-03	∅35×h250 mm	33,300
V03-132-220	4510-04	∅38×h250 mm	42,800
V03-132-226	4510-05	∅45×h350 mm	50,100

### 90° 곡형

Cat. No.	Model	Size	EA Price ₩
V03-132-232	4512-01	∅25×h200 mm	24,800
V03-132-238	4512-02	∅32×h250 mm	28,700
V03-132-244	4512-03	∅35×h250 mm	33,300
V03-132-250	4512-04	∅38×h250 mm	42,800
V03-132-256	4512-05	∅45×h350 mm	50,100



## Large Capacity Vacuum Trap 대용량 진공 트랩, LukeGL® ①

Cat. No.	Model	Flask	Body Size	EA Price ₩
V08-132-433	4529-01	500 ml	∅105×h310 mm	101,200
V08-132-440	4529-02	1000 ml	∅131×h340 mm	134,900

- Borosilicate Glass 3.3
- Trap Stopper Arm : 1/2" OD
- LukeGL®

## Dewer Type Vacuum Trap 드와형 진공 트랩, LukeGL® ②

Cat. No.	Model	Flask	Inner Chamber	EA Price ₩
V08-132-333	4518-11	100 ml	∅40×h60 mm	143,400
V08-132-340	4517-11	500 ml	∅75×h275 mm	286,700
V08-132-347	4517-02	1000 ml	∅85×h300 mm	320,500

- Borosilicate Glass 3.3
- PTFE Stopcock 기본 장착
- LukeGL®

