

## O-Ring Type Vacuum Trap 오링형 진공 트랩

Cat. No.	Model	Body (od×h, mm)	O-ring	Inner Tube (od, mm)	EA Price ₩
V03-132-917	4513-10	36×225	30	12	384,100
V03-132-919	4513-11	46×250	40	16	420,700
V03-132-921	4513-12	56×275	50	25	491,100

- O-ring Joint형으로 Stopper에 Trap Body 연결
- Vacuum Tight, Grease Free, 세척을 위해 쉽게 Stopper와 분리 가능
- Hose Connection : 13.5 mm at Smallest Serration,  
15.5 mm at the largest Serration
- 구성 : Lower Body, Stopper, Viton O-ring, S.S Pinch Clamp
- Chemglass

### Replacement Part

Cat. No.	Model	Description	EA Price ₩
V03-132-951	4513-20	Lower Body, #30 O-Ring	91,600
V03-132-953	4513-21	Lower Body, #40 O-Ring	107,600
V03-132-955	4513-22	Lower Body, #50 O-Ring	137,300
V03-132-957	4513-30	Stopper, #30 O-Ring	181,200
V03-132-959	4513-31	Stopper, #40 O-Ring	190,200
V03-132-961	4513-32	Stopper, #50 O-Ring	231,300
V08-132-626	150-08	Pinch Clamp, #30 O-Ring용	73,600
V08-132-624	150-09	Pinch Clamp, #40 O-Ring용	102,500
V08-132-620	150-21	Pinch Clamp, #50 O-Ring용	136,300
V08-132-912	305-223	Viton O-Ring, #30, (3pk)	18,600
V08-132-913	305-226	Viton O-Ring, #40, (3pk)	21,700
P05-132-579	305-229	Viton O-Ring, #50, (3pk)	25,900



## Schlenk O-ring Type Vacuum Trap 오링형 쉬링크 진공 트랩

Cat. No.	Model	Body (od×h, mm)	O-ring	EA Price ₩
V03-132-901	AF-0200-01	57×300	# 15	737,100

- Borosilicate Glass 3.3
- O-ring Joint형으로 Stopper에 Trap Body 연결
- Two-piece Vacuum Tight, Grease Free, 57mm Flat O-ring Flange
- 구성 : Trap Stopper and Body, #331 Viton o-ring 1개,  
#116 Viton o-ring 2개, Horseshoe Clamp
- Chemglass

### Replacement Part

Cat. No.	Model	Description	EA Price ₩
V08-132-062	AF-0200-02	Trap Body, #15 O-ring	160,500
V03-132-963	AF-0200-03	Stopper, #15 O-ring	383,500
V03-132-916	CG-139-02	57 mm Horseshoe Clamp	116,700
V03-132-644	CG-305-116	#116 Viton O-ring	3,300
V03-132-914	CG-305-331	#331 Viton O-ring	10,100



## Bubbler, Pressure Release 가스 치환 밸브

Cat. No.	Model	Description	EA Price ₩
V03-132-725	4536-01	Rapid Purge Valve, Complete	403,300

- AIRFREE® 또는 Vacuum Aparatus Gas Pursing용
- The check valve prevents oil and air from being pulled into the apparatus
- Hose Connection : od10 mm
- Chemglass

### Replacement Part

Cat. No.	Model	Description	EA Price ₩
V03-132-760	999-02	Spring Loaded PTFE Plug	70,000

