

Vacuum Check Valve 진공 체크 밸브

Cat. No.	Model	Size	Connection	EA Price ₩
V03-132-730	4539-01	∅10×L125 mm	-	48,300
V03-132-740	4540-01	∅10×L100 mm	9 mm Hose Nipple	114,600

- Hose Connection이 있는 타입과 없는 타입 2가지
- Chemglass



Vacuum Peptide Synthesis Vessel 진공 펩티드 합성 반응조

- T-Bore PTFE Stopcock
- Hose connection is sealed on the inner joint for connection to a vacuum source
- MEDIUM or COARSE porosity fritted glass resin support 와 사용 가능
- Stopcock and Hose Connection : od 8 mm
- Chemglass

MEDIUM porosity

Cat. No.	Model	용량	외경	전체높이	Thread	EA Price ₩
V03-132-801	CG-1866-01	10 ml	18 mm	190 mm	GL 14	383,400
V03-132-803	CG-1866-02	25 ml	25 mm	209 mm	GL 25	418,300
V03-132-805	CG-1866-03	50 ml	30 mm	238 mm		453,900
V03-132-807	CG-1866-04	100 ml	51 mm	227 mm	GL 32	492,700
V03-132-809	CG-1866-05	250 ml	57 mm	276 mm		524,700
V03-132-811	CG-1866-06	500 ml	80 mm	300 mm		614,900

COARSE porosity

Cat. No.	Model	용량	외경	전체높이	Thread	EA Price ₩
V03-132-813	CG-1866-10	10 ml	18 mm	190 mm	GL 14	383,400
V03-132-815	CG-1866-11	25 ml	25 mm	209 mm	GL 25	418,300
V03-132-817	CG-1866-12	50 ml	30 mm	238 mm		453,900
V03-132-819	CG-1866-13	100 ml	51 mm	227 mm	GL 32	492,700
V03-132-821	CG-1866-14	250 ml	57 mm	276 mm		524,700
V03-132-823	CG-1866-15	500 ml	80 mm	300 mm		590,700

Replacement Part

Cat. No.	Model	Description	EA Price ₩
B08-132-205	CG-196-01	Solid Cap, PTFE Seal, GL 14 Thread	7,300
B08-132-208	CG-196-03	Solid Cap, PTFE Seal, GL 25 Thread	8,600
B08-132-211	CG-196-04	Solid Cap, PTFE Seal, GL 32 Thread	14,100
S11-132-203	CG-436-67	2mm T-Bore PTFE Plug Only, 15.2/30	47,000

ASTM F128-94 규격에 따른 Porosity

Porosity Grade	ASTM Grade	Pore Size
Medium	M	10 ~ 16 um
Coarse	C	40 ~ 60 um
M	Filtration and Extraction	
C	Gas Dispersion, Gas Washing and Absorption. Coarse Filtration	

