

LC Certified 2ml Vial Kit LC용 2 ml 바이알 키트, 테스트 성적서 포함

Cat. No.	Model	색상	규격	Septa	PK Price ₩
V06-214-301	CERT95-12	Clear	기본형	PTFE/ Silicon	100ea 124,900
V06-214-305	CERT95-12A	Amber			100ea 126,600
V06-214-309	CERT95E-12	Clear	Graduated Marking Spot	PTFE/ Silicon w/Slit	100ea 128,700
V06-214-313	CERT95E-12A	Amber			100ea 132,300
V06-214-317	CERT97-12	Clear	기본형	PTFE/ Silicon w/Slit	100ea 132,300
V06-214-321	CERT97-12A	Amber			100ea 136,700
V06-214-325	CERT97E-12	Clear	Graduated Marking Spot	PTFE/ Silicon w/Slit	100ea 139,200
V06-214-329	CERT97E-12A	Amber			100ea 142,400

- Test Parameter에 대한 성적서가 포함된 바이알과 셉타 일체형 캡 셋트
- Size : Ø12×h32 mm, Neck od9 mm (id6 mm), R.A.M[®]. Vial
- Bonded Royal Blue Open Top Cap with Septa
- 재질
Clear Glass : Type I , Class A, 33 Expansion Clear
Amber Glass : Borosilicate 51A
- Test Parameter
 - Dimensions : Height, Diameter, Bottom Thickness, Neck Length, Thread Profile, Annealing
 - Cap/Septum : Diameter, Height, Thread Integrity, Center Flash, Septa Insertion
 - Chemical Test : HPLC Test-Non-volatile residue, GC Test-Volatile Residue
- JGF (J. G. Finneran)



Mass Spec 2ml Vial Test Kit 질량 분석기용 2 ml 바이알 키트, 테스트 성적서 포함

Cat. No.	Model	색상	규격	Septa	PK Price ₩
V06-214-351	MSC95-12	Clear	기본형	PTFE/ Silicon	100ea 191,300
V06-214-355	MSC95-12A	Amber			100ea 194,800
V06-214-359	MSC95E-12	Clear	Graduated Marking Spot	PTFE/ Silicon w/Slit	100ea 197,200
V06-214-363	MSC95E-12A	Amber			100ea 202,400
V06-214-367	MSC97-12	Clear	기본형	PTFE/ Silicon w/Slit	100ea 207,000
V06-214-371	MSC97-12A	Amber			100ea 210,100
V06-214-375	MSC97E-12	Clear	Graduated Marking Spot	PTFE/ Silicon w/Slit	100ea 212,200
V06-214-379	MSC97E-12A	Amber			100ea 218,800

- Test Parameter에 대한 성적서 포함된 바이알과 셉타 일체형 캡 셋트
- Size : Ø12×h32 mm, Neck od9 mm (id6 mm), R.A.M[®]. Vial
- Bonded Royal Blue Open Top Cap with Septa
- 재질
Clear Glass : Type I , Class A, 33 Expansion Clear
Amber Glass : Borosilicate 51A
- Test Parameter
 - Dimensions : Height, Diameter, Bottom Thickness, Neck Length, Thread Profile, Annealing
 - Cap/Septum : Diameter, Height, Thread Integrity, Center Flash, Septa Insertion
 - GL/MS Background Test
 - LC/MS Background Test
 - Particle Counts : Particle Size, Particle (ml/vial)
- JGF (J. G. Finneran)

