

8-425 Preassembled Cap & Septa 8-425 스크류 캡과 셉타 일체형

• JGF (J. G. Finneran) / Certified (품질 보증서 포함)

with Black PP Open-top Cap

Cat. No.	Model	재질 및 규격	PK Price ₩	CS Price ₩
V04-214-552	804040-08	PTFE / Red Rubber	100ea 33,000	10pk 307,800
V04-214-556	806550-08	Red PTFE / White Silicon	100ea 41,100	10pk 384,200
V04-214-560	804060-08	Red PTFE / White Silicon / Red PTFE	100ea 54,400	10pk 508,200
V04-214-564	806070-08	Red PTFE / White Silicon w / Slit	100ea 52,500	10pk 490,100

with White PP Open-top Cap

Cat. No.	Model	재질 및 규격	PK Price ₩	CS Price ₩
V04-214-566	805050-08W	PTFE / Red Silicon	100ea 43,300	10pk 403,900

Top Hat Cap

Cat. No.	Model	색상	재질 및 규격	PK Price ₩	CS Price ₩
V01-163-196	806550TH-08	Black	Red PTFE /	100ea 55,200	10pk 515,700
V01-163-197	806550TH-08Y	Yellow	White Silicon	100ea 57,900	10pk 541,100

• Top Hat™ : Injection 동안 셉타가 움직이는 것을 방지 하기 위하여 Molded Centery Ring을 가지고 있는 것



8-425 Screw Cap 8-425 스크류 캡

• 재질 : Polypropylene (PP)
• JGF (J. G. Finneran) / Certified (품질 보증서 포함)

Open-top Cap

Cat. No.	Model	색상 / 규격	PK Price ₩	CS Price ₩
V04-214-328	5310L-08	Black	100ea 9,200	10pk 85,200
V04-214-332	5310-08	Black / Shimadzu 전용	100ea 7,000	10pk 64,600
V04-214-336	5310-08B	Blue	100ea 8,100	10pk 75,200
V04-214-340	5310-08G	Green	100ea 8,100	10pk 75,200
V04-214-344	5310-08R	Red	100ea 8,100	10pk 75,200
V04-214-348	5310-08Y	Yellow	100ea 8,100	10pk 75,200
V04-214-352	5310-08W	White	100ea 8,100	10pk 75,200

Closed-top Cap

Cat. No.	Model	색상	규격	PK Price ₩	CS Price ₩
V04-214-458	5360-08	Black	Teflon lined	100ea 32,600	10pk 304,000

• PTFE/F217 Foam Lined



Septa for 8-425 Screw Cap 셉타, 8-425 스크류 캡용

Cat. No.	Model	재질 및 규격	PK Price ₩	CS Price ₩
V01-214-724	601010-08	White PTFE	100ea 5,600	10pk 52,400
V01-214-728	604040-08	PTFE / Red Rubber, 0.045"	100ea 20,500	10pk 191,300
V01-214-730	604550-08	Red PTFE / White Silicon, 0.045"	100ea 22,600	10pk 210,500
V01-214-732	606550-08	Red PTFE / White Silicon, 0.065"	100ea 22,700	10pk 211,700
V01-214-736	606550B-08	Blue PTFE / White Silicon, 0.065"	100ea 31,800	10pk 297,200
V01-214-740	606550T-08	Tan PTFE / White Silicon, 0.065"	100ea 21,200	10pk 197,700
V01-214-744	604060-08	Red PTFE / Silicon / Red PTFE, 0.040"	100ea 38,000	10pk 354,600
V01-214-752	610050-08	PTFE / Silicon, 0.1" Thick	100ea 41,140	10pk 384,450
V01-214-756	606070-08	Red PTFE / White Silicon, w / Slit, 0.060"	100ea 56,400	10pk 526,600

• JGF (J. G. Finneran) / Certified (품질 보증서 포함)



Red PTFE / White Silicon