

## 11 mm Aluminum Seal with Septa 11 mm 알루미늄 씬과 셉타, 2 ml Crimp-top Vial 용

• JGF (J. G. Finneran) / Certified (품질 보증서 포함)

### Silver Aluminum Seal

Cat. No.	Model	Septa 재질	PK Price ₩	CS Price ₩
V01-214-603	5140-11	Clear PTFE / Red Rubber	100ea 19,300	10pk 179,600
V01-214-606	5150-11	Red PTFE / White Silicon	100ea 39,700	10pk 371,000
V01-214-609	5160-11	Red PTFE / White Silicon / Red PTFE	100ea 45,200	10pk 421,600
V01-214-612	5170-11	Red PTFE / White Silicon w/ Slit	100ea 65,400	10pk 610,500
V01-214-618	5110-11	10mil PTFE	100ea 19,600	10pk 183,000
V01-214-619	5150UL-11	1 mm PTFE / Blue Silicon, Ultra Low Bleed	100ea 47,700	10pk 445,400

• Ultra Low Bleed Seal : 5150UL-11

- Very Sensitive Application & Multiple Injection용
- Eliminate Particulates with no Interfering Peaks

### Blue Aluminum Seal

Cat. No.	Model	Septa 재질	PK Price ₩	CS Price ₩
V01-214-621	5140-11B	Clear PTFE / Red Rubber	100ea 23,800	10pk 222,300
V01-214-624	5150-11B	Red PTFE / White Silicon	100ea 45,600	10pk 425,900
V01-214-627	5160-11B	Red PTFE / White Silicon / Red PTFE	100ea 50,100	10pk 468,100

### Green Aluminum Seal

Cat. No.	Model	Septa 재질	PK Price ₩	CS Price ₩
V01-214-630	5140-11G	Clear PTFE / Red Rubber	100ea 23,800	10pk 222,300
V01-214-633	5150-11G	Red PTFE / White Silicon	100ea 45,600	10pk 425,900
V01-214-636	5160-11G	Red PTFE / White Silicon / Red PTFE	100ea 50,100	10pk 468,100

### Red Aluminum Seal

Cat. No.	Model	Septa 재질	PK Price ₩	CS Price ₩
V01-214-639	5140-11R	Clear PTFE / Red Rubber	100ea 23,800	10pk 222,300
V01-214-642	5150-11R	Red PTFE / White Silicon	100ea 45,600	10pk 425,900
V01-214-645	5160-11R	Red PTFE / White Silicon / Red PTFE	100ea 50,100	10pk 468,100



## 11 mm Poly Crimp Seal with Septa 11 mm 플라스틱 씬과 셉타, 2 ml Crimp-top Vial 용

- 모든 11 mm Crimp-top Vial과 같이 사용 가능
- Snap Seal Vial과도 같이 사용 가능
- 기존의 Crimping 또는 Decapping Tool과 같이 사용
- JGF (J. G. Finneran) / Certified (품질 보증서 포함)

### Clear Poly Crimp Seal / 투명 기본형

Cat. No.	Model	Septa 재질	PK Price ₩	CS Price ₩
V01-214-648	5500-11	Molded Septa	100ea 8,600	10pk 80,000
V01-214-651	5540-11	PTFE / Red Rubber	100ea 27,200	10pk 253,400
V01-214-654	5550-11	PTFE / White Silicon	100ea 34,700	10pk 324,300
V01-214-657	5560-11	Red PTFE / White Silicon / Red PTFE	100ea 43,000	10pk 401,100
V01-214-660	5570-11	Red PTFE / White Silicon w/ Slit	100ea 69,800	10pk 651,600
V01-214-663	5500SB-11	Molded Starburst	100ea 8,600	10pk 80,000
V01-214-666	5581-11	Red PTFE / White Silicon w/ Starburst	100ea 68,900	10pk 643,900
V01-214-669	5510-11	10 mil PTFE	100ea 24,000	10pk 224,200

### Red Poly Crimp Seal / 투명 기본형

Cat. No.	Model	Septa 재질	PK Price ₩	CS Price ₩
V01-214-672	5540-11R	PTFE / Red Rubber	100ea 28,700	10pk 268,200
V01-214-675	5550-11R	Red PTFE / White Silicon	100ea 35,600	10pk 332,400
V01-214-678	5560-11R	PTFE / Silicon / PTFE	100ea 44,000	10pk 411,200

