

11 mm Snap Top Cap with Septa 11 mm 스냅-탑 캡과 셉타

- Snap Top Cap 재질 : PE
- Closed Top Cap 중 Starburst가 있는 것은 Needle Injection을 가능하게 함
- Crimper 또는 Decapper의 사용 없이 손으로 작업
- JGF (J. G. Finneran) / Certified (품질 보증서 포함)

Clear Snap Top Cap

Cat. No.	Model	Septa 재질	PK Price ₩	CS Price ₩
V04-214-581	5200-11	w/Molded Septum	100ea 9,600	10pk 89,000
V04-214-582	520070-11	w/Molded Septum, Pre-slit	100ea 21,800	10pk 203,400
V04-214-584	5210-11	10 mil PTFE	100ea 24,000	10pk 224,200
V04-214-588	5200SB-11	Starburst (Closed)	100ea 9,600	10pk 89,400
V04-214-592	5240-11	PTFE / Red Rubber	100ea 35,000	10pk 326,300
V04-214-596	5250-11	PTFE / Silicon	100ea 34,700	10pk 324,300
V04-214-600	5260-11	PTFE / Silicon / PTFE	100ea 42,600	10pk 397,600
V04-214-604	5270-11	PTFE / Silicon w/ Slit	100ea 67,000	10pk 625,800
V04-214-606	5271-11	PTFE / Silicone w/Y-cut	100ea 66,330	10pk 619,850
V04-214-608	5275-11	PTFE / Silicone w/Cross Slit	100ea 66,330	10pk 619,850
V04-214-609	5280-11	PTFE / Silicone/ PTFE w/ Starburst	100ea 88,900	10pk 830,000
V04-214-610	5281-11	PTFE / Silicone w/Starburst	100ea 68,800	10pk 643,000
V04-214-611	5282-11	10mil PTFE w/Starburst	100ea 54,300	10pk 507,000
V04-214-612	5290-11	Viton Lined	100ea 103,900	10pk 970,700



5200SB-11

Blue Snap Top Cap

Cat. No.	Model	Septa 재질	PK Price ₩	CS Price ₩
V04-214-616	5240-11B	PTFE / Red Rubber	100ea 35,000	10pk 326,300
V04-214-620	5250-11B	Red PTFE / White Silicon	100ea 35,600	10pk 332,400
V04-214-624	5260-11B	Red PTFE / White Silicon / Red PTFE	100ea 43,400	10pk 405,300

11 mm GC Snap Top Cap with Septa 11 mm GC용 스냅-탑 캡과 셉타

Cat. No.	Model	색상	Septa 재질	PK Price ₩	CS Price ₩
V04-214-652	5210-11PP	Clear	10mil PTFE	100ea 24,600	10pk 229,400
V04-214-656	5240-11PP		PTFE / Red Rubber	100ea 35,200	10pk 328,200
V04-214-660	5250-11PP		PTFE / Silicon	100ea 39,500	10pk 369,000
V04-214-664	5260-11PP		PTFE / Silicon / PTFE	100ea 47,800	10pk 446,100
V04-214-666	5290-11PP		Black Viton	100ea 104,900	10pk 979,900

- Snap Top Cap 재질 : PP
- Volatile Sample 용으로 Liquid Airtight Seal 제공
- JGF (J. G. Finneran) / Certified (품질 보증서 포함)



Injection Port Septa 인젝션 포트 셉타, 11 mm Crimp-top용

Cat. No.	Model	색상	사용온도	PK Price ₩
V06-214-264	606550BTO-11	Red	400°C	25ea 50,100
V06-214-268	606550GR-11	Green	300°C	50ea 85,200
V06-214-272	606550Y-11	Yellow	350°C	25ea 61,100
V06-214-276	606550LB-11G	Gray	200°C	24ea 22,300
V06-214-280	606550LB-11R	Red	300°C	25ea 40,800
V06-214-282	606MSLB-11	Yellow	250°C	25ea 32,000

- 11 mm Crimp-top용
- JGF (J. G. Finneran)



Model	Specification
606550BTO-11	GC & GC/M Applications, Lowest Bleed Performance / Long Life - Up to 450 Injections
606550GR-11	Low Bleed Performance, Long Life - Up to 100 Injections
606550Y-11	Very Low Bleed Performance, Long Life - Up to 200 Injections
606550LB-11G	Low Bleed Performance, Good Resealing Capability / Long Life - Up to 100 Injections
606550LB-11R	Low Bleed Performance, Good Resealing Capability / Long Life - Up to 100 Injections
606MSLB-11	Low Bleed Yellow PTFE / Cream Silicon Septa