



### 13-425 Screw Cap for 4 ml Vial 13-425 스크류 캡, 4 ml 바이알 용

- 재질 : PP (Polypropylene)
- JGF (J. G. Finneran) / Certified (품질 보증서 포함)

#### Open-top Cap

Cat. No.	Model	색상	PK Price ₩	CS Price ₩
V04-214-380	5310-13	Black	100ea 7,900	10pk 73,000
V04-214-384	5310-13B	Blue	100ea 8,600	10pk 79,600
V04-214-388	5310-13G	Green	100ea 8,600	10pk 79,600
V04-214-392	5310-13R	Red	100ea 8,600	10pk 79,600
V04-214-396	5310-13W	White	100ea 8,600	10pk 79,600
V04-214-400	5310-13Y	Yellow	100ea 8,600	10pk 79,600

#### Closed-top Cap

Cat. No.	Model	색상	규격	PK Price ₩	CS Price ₩
V04-214-466	5360-13	Black	Teflon Lined	100ea 34,400	10pk 321,200

#### Clear Cap, w/10mil Molded Septa

Cat. No.	Model	색상	규격	PK Price ₩	CS Price ₩
V04-214-270	5330-15	Clear	15-425mm	100ea 15,900	10pk 148,600

- 바늘 침투가 용이한 Septa



5360-13



5330-15

### Septa for 13-425 Screw Cap 셉타, 13-425 스크류 캡용

Cat. No.	Model	Septa 재질	PK Price ₩	CS Price ₩
V01-214-784	601010-13	10 mil White PTFE	100ea 6,800	10pk 62,800
V01-214-788	606040-13	PTFE / Red Rubber	100ea 26,300	10pk 245,400
V01-214-792	606050-13	Red PTFE / White Silicon	100ea 33,300	10pk 311,200
V01-214-796	606050G-13	Gray PTFE / White Silicon	100ea 36,400	10pk 339,900
V01-214-800	604060-13	Red PTFE / Silicon / Red PTFE	100ea 45,600	10pk 425,900
V01-214-804	606070-13	Red PTFE / White Silicon W/ slit	100ea 65,100	10pk 607,700

- JGF (J. G. Finneran) / Certified (품질 보증서 포함)



Red PTFE / White Silicon

### 13-425 Preassembled Cap and Septa 13-425 스크류 캡과 셉타 일체형

- JGF (J. G. Finneran) / Certified (품질 보증서 포함)

#### Black PP Open-top Cap

Cat. No.	Model	Septa	PK Price ₩	CS Price ₩
V04-214-668	804060-13	Red PTFE / Silicon / Red PTFE	100ea 64,600	10pk 603,400
V04-214-672	806550-13	Red PTFE / White Silicon	100ea 52,600	10pk 491,300

#### Black Top Hat Cap

Cat. No.	Model	Septa	PK Price ₩	CS Price ₩
V04-214-676	806550TH-13	PTFE / Silicon	100ea 68,300	10pk 637,500

- Top Hat™ : Injection 동안 셉타가 움직이는 것을 방지 하기 위하여 Molded Centery Ring을 가지고 있는 것

