

4 ml Crimp Top & Snap Ring Vial 4 ml 크림프-탑 과 스냅 캡 겸용 바이알, 13 mm Crimp-top

Cat. No.	Model	색상	규격	PK Price ₩	CS Price ₩
V01-214-412	33813-1545	Clear	기본형	100ea 28,800	10pk 269,000
V01-214-424	33813-1545Z		Silanized	100ea 84,400	10pk 788,400
V01-214-428	33813-1545A	Amber	기본형	100ea 33,500	10pk 312,400
V01-214-436	33813E-1545A		눈금 & Marking Spot	100ea 38,700	10pk 361,100
V01-214-440	33813-1545AZ		Silanized	100ea 90,310	10pk 843,150

- Waters WISP 48 position, 717, Breeze와 Shimadzu Autosampler에 적합
- Step Vial Design으로 Limited Volume Insert를 정확하게 Center에 위치 시킴
- Size : Ø15×h45 mm, Ø13 mm Aluminum Seal 및 Snap Cap용
- Silanized Vial : 비활성 처리 (Deactivated Treatment)
- 재질 - Clear Glass : USP Type 1 Borosilicate Glass 3.3
- Amber Glass : USP Type 1 Borosilicate Glass 5.1
- JGF (J. G. Finneran) / Certified (품질 보증서 포함)



- Glass Limited Volume Insert와 Plastic Limited Volume Insert (PP)는 4 ml Screw-top Vial 참조 (1228 page)



13 mm Aluminum Seal with Septa 13 mm 알루미늄 씬과 셉타

Cat. No.	Model	Septa	PK Price ₩	CS Price ₩
V01-214-684	5140-13	PTFE / Red Rubber	100ea 45,300	10pk 422,700
V01-214-687	5150-13	Red PTFE / White Silicon	100ea 52,900	10pk 494,300
V01-214-690	5160-13	Red PTFE / White Silicon / Red PTFE	100ea 60,200	10pk 562,700
V01-214-693	5110-13	10 mil PTFE	100ea 30,700	10pk 286,900

- JGF (J. G. Finneran) / Certified (품질 보증서 포함)



13 mm Snap Top Cap 13 mm 스냅 탑 캡

- JGF (J. G. Finneran) / Certified (품질 보증서 포함)

Clear Snap Top Cap with Septa ①

Cat. No.	Model	Septa	PK Price ₩	CS Price ₩
V04-214-680	5240-13	PTFE / Red Rubber	100ea 43,000	10pk 401,800
V04-214-692	5210-13	10 mil PTFE	100ea 32,200	10pk 300,200



Clear Snap Top Cap with Starburst ②

Cat. No.	Model	Septa	PK Price ₩	CS Price ₩
V04-214-697	5200SB-13	13 mm, with Starburst	100ea 7,800	10pk 73,700

- 캡 상부에 별모양의 틈이 있어 Needle Injection 및 Resealability를 향상 시킴

