

1 ml Shell Vial and Plug Cap, 8 mm Plug 1 ml 셸 바이알과 캡, for Waters WISP®

- Waters WISP™96 Position 및 Shimadzu Autosampler에 적합
- Vial과 Cap 분리 포장
- Type Borosilicate 3.3 or 5.1 Glass 및 PP
- JGF (J. G. Finneran) / Certified (품질 보증서 포함)

1ml Shell Vial, 8×40 mm

1ml Shell Vial with 8 mm Starburst Plug Cap

Cat. No.	Model	재질	색상	PK Price ₩	CS Price ₩
V01-214-454	9800-840	Glass	Clear ①	100ea 22,100	10pk 206,600
V01-214-458	9800-840A		Amber	100ea 31,800	10pk 297,200
V01-214-462	9800P-840	PP	Clear	100ea 30,600	10pk 285,400

0.7 ml Limited Volume PP Shell Vial with Starburst Plug Cap ②

Cat. No.	Model	재질	색상	PK Price ₩	CS Price ₩
V01-214-466	9807P-840	PP	Clear	100ea 35,600	10pk 332,400

1 ml Shell Vial Only

Cat. No.	Model	재질	색상	PK Price ₩	CS Price ₩
V01-214-484	4100-840	Glass	Clear	100ea 16,300	10pk 151,800
V01-214-488	4100-840A		Amber	100ea 18,600	10pk 173,300
V01-214-492	4100P-840	PP	Clear	100ea 18,100	10pk 169,000

0.7 ml PP Shell Vial Only

Cat. No.	Model	재질	색상	PK Price ₩	CS Price ₩
V01-214-494	4107P-840	PP	Clear	100ea 27,000	10pk 252,200
V04-214-150	4107P-840A		Amber ③	100ea 28,400	10pk 264,700

0.75 ml Shell Vial, 8×30 mm

0.75 ml Shell Vial with Starburst Plug Cap

Cat. No.	Model	재질	색상	PK Price ₩	CS Price ₩
V01-214-471	9800-830	Glass	Clear	100ea 26,500	10pk 247,000

Shell Vial Only, 8×30 mm

Cat. No.	Model	용량	재질	색상	PK Price ₩	CS Price ₩
V01-214-473	4100-830	750 ul	Glass	Clear	100ea 17,900	10pk 166,600
V01-214-475	4100P-830	900 ul	PP	Natural	100ea 15,600	-

8mm Snap Plug Cap Only, Starburst Type

Cat. No.	Model	색상	규격	PK Price ₩	CS Price ₩
V01-214-924	5400SB-08	Natural	8 mm	100ea 6,700	10pk 62,400
V01-214-928	5400SB-08B	Blue		100ea 7,100	10pk 65,900
V01-214-932	5400SB-08R	Red		100ea 7,100	10pk 65,900
V01-214-936	5400SB-08Y	Yellow		100ea 7,100	10pk 65,900
V01-214-938	5400TP-08	Natural	8 mm, Easy Pierce	100ea 7,900	10pk 73,400
V01-214-940	540030-08	Natural	8 mm, Silicon Lined	100ea 46,600	10pk 435,500

- 캡상부에 별모양의 틈이 있어 Needle Injection 및 Resealability를 향상시킴
- 재질 : PE (Polyethylene)
- Silicon Lined Plug Cap은 Needle Penetration후 기본형 보다 나은 Resealability를 가짐

Glass limited Volume Insert

Cat. No.	Model	용량	규격	PK Price ₩
V04-214-910	4015BS-529	150 ul	Bottom Spring	100ea 76,700
V04-214-964	402BS-536	200 ul ④		100ea 83,400
V04-214-812	4025-629	250 ul	No Spring ⑤	100ea 104,700
V04-214-816	4025BS-629		Bottom Spring ⑥	100ea 112,500

- Spring 재질 : PE (Polyethylene)

