

2 ml Shell Vial and Plug Cap 2 ml 셸 바이알과 캡, 12×32 mm

- Type I Borosilicate 3.3 및 5.1 Glass 및 PP
- Vial과 Cap 분리 포장
- JGF (J. G. Finneran) / Certified (품질 보증서 포함)

2 ml Shell Vial with 12 mm Conical Snap Plug Cap

Cat. No.	Model	재질	색상	PK Price ₩	CS Price ₩
V05-214-004	9800-1232	Glass	Clear	100ea 35,600	10pk 332,400
V05-214-008	9800-1232A		Amber	100ea 42,500	10pk 396,700
V05-214-012	9800P-1232	PP	Clear	100ea 38,900	10pk 363,500

12 mm Plug Cap Only, Starburst Type

Cat. No.	Model	색상	규격	PK Price ₩	CS Price ₩
V06-214-504	5405SB-12	Natural	12 mm	100ea 14,700	10pk 136,700

- 캡상부가 별 모양의 틈이 있어 Needle Injection 및 Resealability를 향상시킴

250 ul Limited Volume Insert

Cat. No.	Model	용량	재질	규격	PK Price ₩	CS Price ₩
V04-214-816	4025BS-629	250 ul	Glass	Bottom	100ea 112,500	-
V04-214-832	4025PBS-631		PP	Spring	100ea 29,600	10pk 276,700



4ml Shell Vial and Plug Cap 4ml 셸 바이알 과 캡, 15×45mm

- Type I Borosilicate Glass 또는 PP
- JGF (J. G. Finneran) / Certified (품질 보증서 포함)

4 ml Shell Vial with 15 mm Starburst Plug Cap

Cat. No.	Model	재질	색상	PK Price ₩	CS Price ₩
V05-214-016	9800-1545	Glass	Clear	100ea 38,900	10pk 363,100
V05-214-020	9800-1545A		Amber	100ea 43,400	10pk 405,300
V05-214-024	9800P-1545	PP	Clear	100ea 41,700	10pk 389,300

3.0 ml Plastic Limited Volume Vial

Cat. No.	Model	재질	색상	PK Price ₩	CS Price ₩
V05-214-504	9830P-1545	PP	Clear	100ea 66,100	10pk 617,500

15 mm Plug Cap Only, Starburst Type

Cat. No.	Model	색상	규격	PK Price ₩	CS Price ₩
V06-214-508	5405SB-15	Natural	15 mm	100ea 15,400	10pk 143,900

- 캡상부에 별 모양의 틈이 있어 Needle Injection 및 Resealability를 향상시킴
- Spring 재질 : PE (Polyethylene)

15mm Polyethylene Snap Plug


Cat. No.	Model	색상	규격	PK Price ₩	CS Price ₩
V04-214-230	5410SB-15	Clear	15mm	100ea 16,300	10pk 151,700

- Pre-scored Starburst Type : 바늘 침투가 좀 더 용이함

300ul Limited Volume Insert, 6×38mm

Cat. No.	Model	용량	재질	규격	PK Price ₩
V04-214-920	403BS-638X	300 ul	Glass	Bottom	100ea 81,700
V04-214-952	403PBS-638		PP	Spring	100ea 28,700

- Self-Centering 기능을 가지고 있음



5410SB-15