

Headspace Vial, 20 mm Crimp Style 20 mm 용 크림프 탑 바이알

- Type I Borosilicate 3.3 Clear Glass
- JGF (J. G. Finneran) / Certified (품질 보증서 포함)

Flat Bottom Headspace Clear Vial, Crimp Top ①

Cat. No.	Model	용량	Size (Ø×h)	PK Price ₩	CS Price ₩
V05-214-182	36020-2238	6 ml	22×38 mm	100ea 64,800	10pk 605,100
V05-214-186	310020-2346	10 ml	23×46 mm	100ea 64,500	10pk 602,200
V05-214-190	320020-2375	20 ml	23×75 mm	100ea 65,500	10pk 612,200
V05-214-194	327020-3060	27 ml	30×60 mm	100ea 151,100	10pk 1,411,900

- Heating 효율을 최대화할 수 있음
- Agilent, Carlo Erba, Perkin Elmer and Shimadzu등과 같이 사용 가능
- 320020E-2375 (Graduation & Marking Spot), 320020M-2375 (Marking Spot)

Flat Bottom Headspace Vial, 18mm Thread ②

Cat. No.	Model	용량	Size (Ø×h)	색상	PK Price ₩	CS Price ₩
V05-214-212	310018F-2346	10 ml	23×46 mm	Clear	100ea 61,300	10pk 572,500
V05-214-216	31018F-2346A			Amber	100ea 88,200	10pk 823,900
V05-214-220	320018F-2375	20 ml	23×75 mm	Clear	100ea 67,500	10pk 630,100
V05-214-224	32018F-2375A			Amber	100ea 90,400	10pk 844,200

- Heating 효율을 최대화할 수 있음
- Agilent, Carlo Erba, Perkin Elmer and Shimadzu등과 같이 사용 가능

Round Bottom Headspace Clear Vial, Crimp Top ③

Cat. No.	Model	용량	Size (Ø×h)	PK Price ₩	CS Price ₩
V05-214-200	310020A-2346	10 ml	23×46 mm	100ea 75,200	10pk 702,200
V05-214-204	320020A-2375	20 ml	23×75 mm	100ea 77,000	10pk 719,500

- 바이알 내부 압력에 더 잘 견딤 (고압용)
- Perkin-Elmer, Tekmar, L.E.A.P. Technology 및 Varian등과 같이 사용 가능

Round Bottom Headspace Vial, 18 mm Thread ④

Cat. No.	Model	용량	Size (Ø×h)	색상	PK Price ₩	CS Price ₩
V05-214-232	310018R-2346	10 ml	23×46 mm	Clear	100ea 75,900	10pk 709,300
V05-214-236	31018R-2346A			Amber	100ea 84,200	10pk 786,100
V05-214-240	320018R-2375	20 ml	23×75 mm	Clear	100ea 80,100	10pk 748,400
V05-214-244	32018R-2375A			Amber	100ea 92,000	10pk 859,800

- Heating 효율을 최대화 할 수 있음
- Agilent, Carlo Erba, Perkin Elmer and Shimadzu등과 같이 사용 가능

SPME Specialty Vial ⑤

Cat. No.	Model	용량	Size (Ø×h)	PK Price ₩	CS Price ₩
V05-214-247	31020SP-2346	10 ml	23×46 mm	100ea 78,200	10pk 730,800
V05-214-249	32020SP-2375	20 ml	23×75 mm	100ea 78,200	10pk 730,800

- 일반 Headspace Vial보다 Thicker Tip을 가지고 있는 바이알
- Thicker Tip (Thin Seal)은 Tight Seal을 유지시킴
- For SPME applications, a 20 mm aluminum seal with a thin septum would be ideal in order extend the life of the needle
- Agilent, Carlo Erba, Perkin Elmer and Shimadzu등과 같이 사용 가능

Metal Threaded Closure, 18 mm Silver ⑥

Cat. No.	Model	Septa	PK Price ₩	CS Price ₩
V05-214-270	5350-18M	0.05" PTFE/Silicon Liner	100ea 136,300	10pk 1,273,000
V05-214-274	5356B-18M	Blue 0.06" PTFE/Silicon Liner	100ea 138,800	10pk 1,296,600
V05-214-278	5351TB-18M	Translucent Blue 0.125" PTFE/Silicon Liner	100ea 142,300	10pk 1,329,200
V05-214-282	5001-20	20 mm Butyl Stopper	100ea 43,300	10pk 404,700



다음 페이지에서 계속...