

Headspace Vial, 20 mm Crimp Style 에서 계속...
20 mm Aluminum Seal with Septa / 셉타 일체형 알루미늄 씰
Open Center Hole Seal with Septa, Silver

Cat. No.	Model	Septa	PK Price ₩	CS Price ₩
V05-214-724	5110-20	10 mil PTFE	100ea 38,600	10pk 360,200
V05-214-728	5140-20	PTFE / Gray Butyl, ①	100ea 122,800	10pk 1,147,700
V05-214-730	5142-20	PTFE / Gray Butyl Flange	100ea 168,200	10pk 1,571,900
V05-214-732	5150-20	PTFE / Silicon, ②	100ea 77,400	10pk 722,900
V05-214-734	5150UL-20	3 mm Thick PTFE / Blue Silicon, Ultra Low Bleed	100ea 104,100	10pk 972,400

- One-piece Unlined Aluminum Seal
- Ultra Low Bleed Septa : 5150UL-20
 - Very Sensitive Application & Multiple Injection용
 - Eliminate Particulates with no Interfering Peaks

Pressure Release Seal with Septa, Silver

Cat. No.	Model	Septa	PK Price ₩	CS Price ₩
V05-214-736	5141PR-20	PTFE / Gray Butyl ③	100ea 198,600	10pk 1,856,100
V05-214-740	5142PR-20	PTFE / Gray Butyl Rubber	100ea 231,600	10pk 2,164,000
V05-214-744	5143PR-20	PTFE / Silicon	100ea 191,000	10pk 1,785,100

- 바이알 내부 압력이 3.0±0.5 bar를 초과하면 씰이 오픈 됨

Tin Plate Magnetic Seal with Septa ④

Cat. No.	Model	Septa	PK Price ₩
V05-214-748	5144TP5-20	0.131" PTFE / Gray Butyl, 5 mm Hole	100ea 187,900
V05-214-750	5144TP8-20	0.131" PTFE / Gray Butyl, 8 mm Hole	100ea 188,200
V05-214-752	5145TP5-20	0.1" PTFE / Silicon, 5 mm Hole	100ea 168,200
V05-214-753	5145TP8-20	0.1" PTFE / Silicon, 8 mm Hole	100ea 171,100

- Magnetic Vial Transport를 이용하는 Autosampler용

Tinner Seal for SPME Application ⑤

Cat. No.	Model	Septa	PK Price ₩
V05-214-771	5145TP8-20SP	1.5 mm (0.06") Blue PTFE / White Silicon	100ea 200,500
V05-214-773	5156TP8-20SP	1.5 mm (0.06") Natural PTFE / Silicon	100ea 110,000

- 20 mm Gold Magnetic Seal with 8 mm Hole

BiMetal Magnetic Seal ⑥

Cat. No.	Model	색상	Septa	PK Price ₩	CS Price ₩
V05-214-755	5150BM00-20R	Red	0.1" PTFE / Silicon	100ea 94,400	10pk 882,200
V05-214-758	5150BM00-20B	Blue	0.1" PTFE / Silicon	100ea 94,400	10pk 882,200
V05-214-762	5150BM18-20R	Red	0.118" PTFE / Silicon	100ea 112,500	10pk 1,050,800
V05-214-766	5150BM18-20B	Blue	0.118" PTFE / Silicon	100ea 112,500	10pk 1,050,800
V05-214-768	5150BMUL-20B	Blue	3 mm Thick PTFE / Blue Silicon, Ultra Low Bleed	100ea 117,300	10pk 1,095,500

- Outer Material : Blue or Red Aluminum • Inner Material : Silver Tin
- Aluminum Seal보다 Crimping을 쉽게 함





20 mm Aluminum Seal Only / 알루미늄 씰

Cat. No.	Model	Type	색상	PK Price ₩	CS Price ₩
V06-214-164	5100-20	Open-Hole	Silver	100ea 10,200	10pk 94,900
V06-214-168	5100-20R		Red	100ea 10,340	10pk 97,240
V06-214-176	5100-20M	Tin Plate	-	100ea 28,270	10pk 263,890
V06-214-188	5100CT-20	Center Tear	Silver	100ea 9,200	10pk 85,900
V06-214-192	5100CT-20R		Red	100ea 9,130	10pk 85,360
V06-214-216	5100TA-20	Total Tear	Silver	100ea 9,000	10pk 83,500
V06-214-220	5100TA-20R		Red	100ea 13,200	10pk 125,400

- Open Hole Seal : 씰 상단 중앙이 뚫려 있음
- Center Tear Seal : 씰 중앙 상단의 Disc를 떼어 낼 수 있음
- Total Tear Seal : 씰 중앙 상단 Disc 및 씰 전체를 떼어 낼 수 있음
- 세부적인 품목은 Vial Accessory란 참조

Septa for 20 mm Aluminum Seal, Ø20mm

Cat. No.	Model	Septa	PK Price ₩	CS Price ₩
V05-214-784	601010-20	10 mil PTFE	100ea 17,400	10pk 161,900
V05-214-788	610040G-20	PTFE / Gray Butyl	100ea 100,000	10pk 934,600
V05-214-790	610040G-20M	PTFE / Gray Butyl Flange	100ea 168,200	10pk 1,571,900
V05-214-792	610050-20S	PTFE / Silicon, 10/90 mils	100ea 64,700	10pk 604,600



Stopper for 20 mm Crimp-top Vial은 1255 Page 참조