

Amber Sample Vial with Cap 갈색 샘플 바이알, 캡 일체형

- 재질 : Type I Borosilicate Glass 5.1
- 광원에 민감한 제품의 Sampling 및 저장용
- 개별 분리형 Cellular 청결 포장된 고품질 바이알
- Qorpak

with Teflon Lined Cap Attached

Cat. No.	Model	용량	Size (Ø × h)	Neck	PK Price ₩	
V05-134-552	C-06786	2 ml	12 × 35 mm	8-425	144ea	318,900
V05-134-556	C-06787	4 ml	15 × 45 mm	13-425	144ea	322,000
V05-134-560	C-06788	8 ml	17 × 60 mm	15-425	144ea	307,000
V05-134-564	812020-1965A	12 ml	19 × 65 mm		100ea	130,400
V05-134-568	C-06789	16 ml	21 × 70 mm	18-400	144ea	445,300
V05-134-572	C-01056	20 ml	27.75 × 57.5 mm	24-400	72ea	227,200
V05-134-576	C-01065	40 ml	27.75 × 95 mm		72ea	257,200
V05-134-580	C-01066	60 ml	27.75 × 140 mm		72ea	284,400

- 유기 화합물과 같이 사용하기에 적합
- 812020-1965A는 PP Black Teflon Lined Cap 임

Polyseal Cone Lined Cap Attached

Cat. No.	Model	용량	Size (Ø × h)	Neck	PK Price ₩	
V05-134-612	C-01061	20 ml	27.75 × 57 mm	24-400	72ea	192,900
V05-134-616	C-04892	40 ml	27.75 × 95 mm		72ea	235,900

- 완벽한 Sealing을 보장하며, 내화학성을 가지고 있음

with White Rubber Lined Cap Attached

Cat. No.	Model	용량	Size (Ø × h)	Neck	PK Price ₩	
V05-134-654	C-04918	20 ml	27.75 × 57 mm	24-400	72ea	210,600
V05-134-656	C-05181	40 ml	27.75 × 95 mm		72ea	243,100

- FDA 승인 된 라이너로 식품 및 제약용에 사용하기에 적합, Autoclavable

with PTFE / Silicon Lined Open-top Cap Attached


Cat. No.	Model	용량	Size (Ø × h)	Neck	PK Price ₩	
V05-134-771	812030-1965A	12 ml	19 × 65 mm	15-425	100ea	155,200
V05-134-775	816030-2170A	16 ml	21 × 70 mm	18-400	100ea	175,500
V05-134-780	C-01058	20 ml	27.75 × 57 mm	24-400	72ea	284,300
V05-134-784	C-01075	40 ml	27.75 × 95 mm		72ea	284,900
V05-134-786	C-07867	60 ml	27.75 × 140 mm		72ea	319,000




Teflon Lined Cap



Polyseal Cone Lined Cap



White Rubber Lined Cap



Compound Wide-Mouth Vial with Cap 단형 광구 바이알, 캡 일체형

- 재질 : Type I Borosilicate Glass 3.3 & 5.1
- PTFE Lined Cap, 고품질 바이알
- Cap Size : 13-425 Thread
- Qorpak

Clear Vial with Teflon Lined Cap

Cat. No.	Model	용량	Size (Ø × h)	Neck	PK Price ₩	
V05-134-258	C-08763	1 ml	14.8 × 19 mm	13-425	144ea	263,300
V05-134-262	C-08764	1.3 ml	14.8 × 22 mm		144ea	266,600
V05-134-266	C-08765	2 ml	14.8 × 26 mm		144ea	269,200

Amber Vial with Teflon Lined Cap

Cat. No.	Model	용량	Size (Ø × h)	Neck	PK Price ₩	
V05-134-762	C-08767	1.3 ml	14.8 × 22 mm	13-425	144ea	284,800
V05-134-766	C-08768	2 ml	14.8 × 26 mm		144ea	284,900

