

Clear Sample Vial with Black Teflon/Rubber Lined Cap, Kimble®

투명 샘플 바이알, 내열 & 내화학 라이너

- Borosilicate Glass 3.3, USP Type I, ASTM E438 Type I, Class A
- 바이알과 캡 모두 Autoclave 가능
- PTFE-faced white rubber liner가 부착된 Closed-top Black Phenolic Cap
- DWK Lifescience / Kimble®

캡 분리 포장형 ①

Cat. No.	Model	용량	Size (Ø×h)	캡 규격	PK Price ₩
V05-137-640	60812B-1232	1.8ml	12×32 mm	8-425	200ea 294,800
V05-137-643	60812B-12	2ml	12×35 mm		200ea 307,300
V05-137-646	60812B-1	4ml	15×45 mm	13-425	200ea 333,100
V05-137-649	60812B-2	8ml	17×60 mm	15-425	200ea 394,700
V05-137-652	60812B-3	12ml	19×65 mm		200ea 443,100
V05-137-655	60812B-4	16ml	21×70 mm	18-400	200ea 512,900
V05-137-658	60812B-6	24ml	23×85 mm	20-400	200ea 546,900



60812B-3

캡 일체형 ②

Cat. No.	Model	용량	Size (Ø×h)	캡 규격	PK Price ₩
V05-137-663	60811B-12	2ml	12×35 mm	8-425	288ea 732,300
V05-137-666	60811B-1	4ml	15×45 mm	13-425	144ea 320,700
V05-137-669	60811B-2	8ml	17×60 mm	15-425	144ea 377,200
V05-137-672	60811B-3	12ml	19×65 mm		144ea 400,000
V05-137-675	60811B-4	16ml	21×70 mm	18-400	144ea 481,900
V05-137-678	60811B-5	20ml	28×57 mm	24-400	72ea 320,800
V05-137-681	60811B-6	24ml	23×85 mm	20-400	144ea 550,100
V05-137-684	60811B-612	25ml	28×70 mm	24-400	72ea 340,500
V05-137-687	60811B-10	40ml	28×95 mm		72ea 413,000



- Partitioned Tray에 포장되어짐

Clear Sample Vial with White Teflon/Rubber Lined Cap, Kimble®

투명 샘플 바이알

Cat. No.	Model	용량	Size (Ø×h)	캡 규격	PK Price ₩
V05-137-693	60940A-2	2ml	12×35 mm	8-425	144ea 280,600
V05-137-696	60940A-4	4ml	15×47 mm	13-425	144ea 201,000
V05-137-699	60940A-8	8ml	17×63 mm	15-425	144ea 244,100
V05-137-702	60940A-12	12ml	19×67 mm		144ea 266,100
V05-137-705	60940A-16	16ml	21×72 mm	18-400	144ea 299,900
V05-137-708	60940A-24	24ml	23×87 mm	20-400	144ea 365,100



- Borosilicate Glass 3.3, USP Type I, ASTM E438 Type I, Class A
- 바이알과 캡 모두 Autoclave 가능
- PTFE-faced White Rubber Liner가 부착된 Closed-top White Urea Cap
- Partitioned Tray에 포장되어짐
- DWK Lifescience / Kimble®

EPA Water Analysis Clear Vial with Open Top Cap EPA Water 분석 투명 바이알, Kimble®

Cat. No.	Model	용량	Size (Ø×h)	캡 규격	PK Price ₩
V05-137-773	60961C-4	20ml	28×60 mm	24-400	72ea 311,700
V05-137-776	60961C-6	30ml	28×73 mm		72ea 371,600
V05-137-782	60961C-912	40ml	28×98 mm		72ea 334,500



- Borosilicate Glass 3.3, USP Type I, ASTM E438 Type I, Class A
- Water 샘플들의 수집 및 저장 용도
- 바이알과 캡 모두 Autoclave 가능
- 오염물 분석을 위한 시험 절차 수립 지침 준수 (Ref. EPA Methods 601, 602, 603, 604)
- PTFE-faced(0.005") silicone rubber(0.120") septa가 부착된 Open-Top White Polypropylene Cap
- DWK Lifescience / Kimble®