

High Recovery Vial 하이 리커버리 바이알


- Open-top Cap과 PTFE/Silicon Septa 포함
- Chemglass

High Recovery Vial with Open Top Cap

Cat. No.	Model	용량	Size (Ø×h)	Cap	PK Price ₩
V05-132-080	4915-01	1.5 ml	12×32 mm	9 mm	100ea 349,800
V05-132-180	4916-30	4ml	15×45 mm	13-425	100ea 416,600
V05-132-184	4916-31	7.3ml	17×60 mm	15-425	100ea 549,400
V05-132-188	4916-32	11.4ml	19×70 mm		100ea 583,500
V05-132-192	4916-33	19.2ml	28×57 mm	22-400	100ea 627,700
V05-132-196	4916-34	23.7ml	28×67 mm		100ea 641,800
V05-132-200	4916-35	33ml	28×95 mm	24-400	100ea 764,900

High Recovery Vial Only

Cat. No.	Model	용량	Size (Ø×h)	Cap	PK Price ₩
V05-132-140	4916-01	3.7ml	15×45 mm	13-425	100ea 344,500
V05-132-144	4916-02	7.3ml	17×60 mm	15-425	100ea 466,700
V05-132-148	4916-03	11.4ml	19×70 mm		100ea 488,300
V05-132-152	4916-04	19.2ml	28×57 mm	22-400	100ea 505,700
V05-132-156	4916-05	23.7ml	28×67 mm		100ea 536,100
V05-132-160	4916-06	33ml	28×95 mm	24-400	100ea 628,100



1.5 ~ 10 ml

20 ~ 60 ml



Mini Reaction Vial / V-Vial V 형 바이알

Cat. No.	Model	용량	Size (Ø×h)	Cap	PK Price ₩
V05-132-528	4913-001	1 ml	13×44 mm	13-425	12ea 238,200
V05-132-535	4913-003	3 ml	20×50 mm	20-400	12ea 259,800
V05-132-542	4913-005	5 ml	20×65 mm	20-400	12ea 267,800

- Open-top Cap과 PTFE/Silicon Septa 포함
- Chemglass



2.5ml Clear Vial, 13-425mm Thread 2.5ml 단형 바이알

Cat. No.	Model	색상	Capacity	Size	PK Price ₩	CS Price ₩
V04-214-800	32513-1528	Clear	2.5ml	15×28 mm	100ea 36,900	10pk 344,200

- Type I Borosilicate, 1st Hydrolytic Class Clear Glass
- JGF (J. G. Finneran)



Straight-Sided Vial Only, Kimble® 내부 일자형 바이알, 캡 별도

Cat. No.	Model	용량	Size (od×h)	GPI Finish	PK Price ₩	CS Price ₩
V05-137-131	60690-12	2 ml	13×39 mm	15-425	250ea 687,800	4pk 2,571,100
V05-137-135	60690-1	4 ml	15×45 mm	18-400	250ea 689,500	4pk 2,577,300
V05-137-139	60690-2	7 ml	18×60 mm	20-400	250ea 792,900	4pk 2,964,000
V05-137-143	60690-4	15 ml	21×70 mm	24-400	250ea 861,900	4pk 3,221,800

- 재질 : 33 Expansion, Low Extractable Borosilicate Glass
- 액체 또는 Powder 샘플의 저장, 추출, 분석, 이동시 용이
- 고속 대량 스크리닝에 적합
- 실험용 자동화 기기에서 고성능을 자랑함
- 내부는 Conical Shape으로 Sample 회수에 용이
- Closure은 별도 구매 (1249 Page 참조)
- DWK Lifescience / Kimble®

