

■ Closed Top Micro Vial, Graduated, Kimble® 눈금형 마이크로 바이알

- 내부가 등고선 형태로 되어 있어 제한된 양의 시료 사용에 적합
- 하단이 좁아지는 형태로 Sample Recovery가 효율적임
- Solid-Top PTFE Faced White Rubber Liner Cap 포함
- DWK Lifescience / Kimble®

Clear Micro Vial, 투명 마이크로 바이알

Cat. No.	Model	용량	Size (Ø×h)	GPI Finish	PK Price ₩
V05-137-147	60702-1	1 ml	14×49 mm	13-425	12ea 314,000
V05-137-151	60702-2	2 ml	21×40 mm	20-400	12ea 331,900
V05-137-155	60702-3	3 ml	21×51 mm		12ea 337,700
V05-137-159	60702-5	5 ml	21×62 mm	24-400	12ea 369,600
V05-137-163	60702-10	10 ml	25×73 mm		12ea 476,400

- 재질 : 33 Expansion, Low Extractable Borosilicate Glass

Amber Micro Vial, 갈색 마이크로 바이알

Cat. No.	Model	용량	Size (Ø×h)	GPI Finish	PK Price ₩
V05-137-167	60707-1	1 ml	14×49 mm	13-425	12ea 546,700
V05-137-171	60707-2	2 ml	21×40 mm	20-400	12ea 581,400
V05-137-175	60707-3	3 ml	21×51 mm		12ea 591,700
V05-137-179	60707-5	5 ml	21×62 mm	24-400	12ea 655,300

- 재질 : 51 Expansion Borosilicate Amber Glass



■ Open Top Micro Vial, Kimble® 오픈탑 마이크로 바이알

- 내부가 등고선 형태로 되어 있어 제한된 양의 시료 사용에 적합
- 하단이 좁아지는 형태로 Sample Recovery가 효율적임
- Open Top Cap with PTFE Faced Silicon Septa 포함
- DWK Lifescience / Kimble®

Clear Open Top Micro Vial, 투명 오픈 탑 마이크로 바이알

- 재질 : 33 Expansion, Low Extractable Borosilicate Glass

without Graduation

Cat. No.	Model	용량	Size (Ø×h)	GPI Finish	PK Price ₩
V05-137-183	60710-110	0.1 ml	12×32 mm	8-425	12ea 323,000
V05-137-187	60710-310	0.3 ml	14×36 mm	13-425	12ea 301,100

with Graduation

Cat. No.	Model	용량	Size (Ø×h)	Cap	PK Price ₩
V05-55-602	60700-1	1 ml	14×49 mm	13-425	12ea 330,600
V05-55-606	60700-2	2 ml	21×40 mm	20-400	12ea 349,400
V05-55-610	60700-3	3 ml	21×51 mm		12ea 355,500
V05-55-614	60700-5	5 ml	21×62 mm	24-400	12ea 389,200
V05-55-618	60700-10	10 ml	25×73 mm		12ea 437,800



Amber Open Top Micro Vial, 갈색 오픈 탑 마이크로 바이알

- 재질 : 51 Expansion Borosilicate Amber Glass

without Graduation

Cat. No.	Model	용량	Size (Ø×h)	GPI Finish	PK Price ₩
V05-137-211	60715-110	0.1 ml	12×32 mm	8-425	12ea 455,700
V05-137-215	60715-310	0.3 ml	14×36 mm	13-425	12ea 495,500

with Graduation

Cat. No.	Model	용량	Size (Ø×h)	Cap	PK Price ₩
V05-55-642	L60705-1	1 ml	14×49 mm	13-425	12ea 546,700
V05-55-646	L60705-2	2 ml	21×40 mm	20-400	12ea 581,400
V05-55-650	L60705-3	3 ml	21×51 mm		12ea 591,700
V05-55-654	L60705-5	5 ml	21×62 mm	24-400	12ea 655,300

