

Crimp Top Micro Vial Only, Kimble® 크림프 탑 마이크로 바이알, 캡 별도

- 재질 : 33 Expansion, Low Extractable Borosilicate Glass
- 내부가 V형태로 되어 있어 주사기 바늘로 남은 잔여 샘플 채취에 용이
- Closure는 별도 구매 (1253 ~ 1255 Page 참조)
- DWK Lifescience / Kimble®

Ungraduated

Cat. No.	Model	용량	Size	A/S Finish	PK Price ₩
V05-137-111	60730-310	0.3 ml	od14×h36 mm	13	12ea 295,500

Graduated

Cat. No.	Model	용량	Size	A/S Finish	PK Price ₩
V05-137-115	60720-1	1 ml	od14×h49 mm	13	12ea 331,900
V05-137-119	60720-2	2 ml	od21×h40 mm	20	12ea 351,400
V05-137-123	60720-3	3 ml	od21×h51 mm		12ea 378,600
V05-137-127	60720-5	5 ml	od21×h62 mm		12ea 416,400



Microwave Reaction Vial Package 마이크로 웨이브용 바이알 팩

- Heavy-Wall, Type I, Class A Borosilicate Glass
- High Resistance to Both Chemical and Thermal Shock
- Max Pressure : up to 300 psi (20 bar)
- Complete Package : Microwave Vial, 20mm Aluminum Seal With 0.125" PTFE/Silicon Septa, PTFE Coated Stirring Bar
- Chemglass

Complete Package

Cat. No.	Model	용량	Size(od x h)	Bottom	PK Price ₩
V05-132-570	4920-01	2 - 5 ml	Ø16×h83 mm	Round	100ea 1,035,400
V05-132-573	4920-02	10 - 20 ml	Ø28×h83 mm		100ea 3,019,900
V05-132-576	4920-03	0.5 - 2 ml	Ø16×h83 mm	Tapered	100ea 1,317,300
V05-132-579	4920-04	0.2 - 0.5 ml	Ø8×h83 mm	-	100ea 2,105,900

- 4920-04 Microwave Vial은 Crimp Seal & Septa와 연결하기 위하여 PEEK Adapter Collar가 필요함 (Not Included)

Spare Accessory

Cat. No.	Model	Description	PK Price ₩
V05-132-600	4920-10	20ml Aluminum Seal with Blue PTFE/Silicon Septa	200ea 462,800
V05-132-605	4920-15	Stir Bar for 4920-01, Ø3×L10mm	50ea 326,700
V05-132-610	4920-16	Stir Bar for 4920-02, Ø6×L15mm, Egg Shape	12ea 62,700
V05-132-615	4920-17	Stir Bar for 4920-03, w5.25×h11mm, Triangular	50ea 375,200
V05-132-620	4920-18	Stir Bar for 4920-04, Ø2×L5mm	50ea 276,600



Rotary Evaporator Adapter for Microwave Vial 회전농축기 바이알 어댑터, 마이크로 웨이브 바이알용

Cat. No.	Model	Joint	Flange	EA Price ₩
V05-132-590	1318-A-10	24/40	20mm Crimp	305,100
V05-132-593	1318-A-12	29/42		318,900

- 구성 : PTFE Joint Adapter, Thread Black Acetal Compression Cap, #017 Viton O-ring
- Chemglass

Spare Viton O-ring

Cat. No.	Model	Description	PK Price ₩
V05-132-625	305-017	Viton O-ring, #017	12ea 19,000



- Microwave Vial용 Heating Block은 591 Page 참조

