


GLS80 Screw Cap with Hose Connection GLS80 호스 연결용 스크류 캡

Cat. No.	Model	연결 포트	Connector 외경	Bore	EA Price ₩
B08-460-110	D804-06	2	od 9 ~ 12 mm	6 mm	149,200
B08-460-112	D804-12	3			171,900

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- GL80 광구병용 Flexible Distributor
- Hose Connector : Bent (Upper), Straight (Lower)
- Upper Hose Connector : 2 Step (od 1 : 9mm, od 2 : 12mm)
- BOLA



GLS80 Screw Cap with Screw Fitting GLS80 스크류 피팅용 스크류 캡

- GL80 광구병용 Multiple Distributor
- 재질 : PP & PTFE, Autoclavable 121°C, 사용 온도 : 0 ~ 110°C
- GL Screw 나사선용 스크류 피팅은 123 Page 참조
- BOLA

PP Cap & PTFE Body

Cat. No.	Model	연결 포트	Neck	EA Price ₩
B08-460-114	D754-08	6	GL14 x 6	299,200
B08-460-116	D754-16	4	GL18 x 4	230,700
B08-460-118	D754-24	4	GL18 x 3 / GL25 x 1	240,100
B08-460-120	D754-36	3	GL18 x 2 / GL32 x 1	376,600

PP Cap & PP Body

Cat. No.	Model	연결 포트	Neck	EA Price ₩
B08-460-122	D750-16	4	GL18 x 4	234,800
B08-460-124	D750-24	4	GL18 x 3 / GL25 x 1	196,500



D754-08




D754-16




D754-24



D754-36



D750-16



D750-24

GLS80 Distributor for Reaction Vessel GLS80 광구병용 반응기 커버

- GLS80 광구병을 반응기로 사용하기 위한 캡
- 재질 : PP & PTFE, Autoclavable 121°C, 사용 온도 : 0 ~ 110°C
- GL Screw 나사선용 스크류 피팅은 123 Page 참조
- BOLA

PP Cap & PTFE Body

Cat. No.	Model	Center	Side	EA Price ₩
B08-460-130	D748-16	29/32	GL18 x 4	542,500
B08-460-132	D748-40	29/32	GL18 x 2 / 29/32 x 2	670,200
B08-460-134	D748-60	29/32	14/23 x 2 / 29/32 x 2	830,800

PP Cap & PP Body

Cat. No.	Model	Center	Side	EA Price ₩
B08-460-136	D746-16	29/32	GL18 x 4	351,900

PP Cap & PTFE Body

Cat. No.	Model	Center / Shaft 지름	Side	EA Price ₩
B08-460-138	D744-16	GL25 / 8mm	GL18 x 4	529,200
B08-460-140	D744-24	GL25 / 10 mm	GL18 x 4	535,500


- Exchangeable Shaft Guide와 고온용 PPS Screw Cap 포함



D748-16



D748-60



D746-16



D744-16



D744-24

