

Open Top Cap with Septa 일체형 오픈탑 캡과 셉타

Cat. No.	Model	Cap 규격	PK Price ₩
V05-216-664	99409	8-425	50ea 70,900
V05-216-666	99419	13-425	50ea 81,510
V05-216-669	99429	15-425	50ea 113,520
V05-216-672	99439	20-400	50ea 133,100
V05-216-678	99459	24-400	50ea 143,220

- Phenolic Open-top Cap과 PTFE/Silicon Septa 일체형
- LK Labkorea



Glass Snap Cap Vial 유리 스냅 바이알

Cat. No.	Model	용량	Size (Ø × h)	PK Price ₩
V06-263-604	60975L-1	4 ml	15 × 45 mm	72ea 47,300
V06-263-612	60975L-3	12 ml	21 × 50 mm	72ea 77,200
V06-263-616	60975L-4	16 ml	25 × 52 mm	72ea 93,500
V06-263-620	60975L-5	20 ml	27 × 55 mm	72ea 109,200

- 재질 : 51 Expansion Borosilicate Glass
- PE 재질의 Closure가 포함되어 있음
- DWK Lifescience / Kimble®



Polypropylene Hinged Vial 뚜껑 일체형 플라스틱 바이알

Cat. No.	Model	용량	Size (Ø × h)	PK Price ₩
V06-134-559	P-10049	75 ml	50 × 50 mm	100ea 178,000
V06-134-563	P-10052	90 ml	33 × 93 mm	100ea 131,200
V06-134-567	P-09773	120 ml	43 × 90 mm (장형)	100ea 81,300
V06-134-571	P-10055		63 × 50 mm (단형)	100ea 201,400
V06-134-575	P-10045	240 ml	75 × 75 mm	100ea 99,400

- 120 ml와 240 ml는 Graduation(눈금)을 가지고 있음
- 요청에 의해 멸균된 제품도 공급 가능
- Qorpak



Storage Vial 스크류캡 바이알

- 재질 : PP Tube, HDPE Closure
- Diagnostic Specimen 또는 Other Solution의 저장 및 이송에 적합
- Tarsons

Sterile, 멸균 바이알

Cat. No.	Model	용량	Size (Ø × h)	PK Price ₩
C25-46-511	T523080	5 ml	14.5 × 64 mm	500ea 89,800
C25-46-513	T523160	10 ml	14.5 × 84 mm	500ea 125,800

Non-Sterile, 비멸균 바이알

Cat. No.	Model	용량	Size (Ø × h)	PK Price ₩
C25-46-521	T523070	5 ml	14.5 × 64 mm	500ea 75,400
C25-46-523	T523150	10 ml	14.5 × 84 mm	500ea 104,200



Cryogenic Vial / 냉동 바이알은 1136 Page 참조

