

## High Purity Teflon Sample Mini Vial 고순도 테프론 샘플 미니 바이알

- 재질 : PFA (Perfluoroalkoxy), 250°C 내열
- 불산 사용이 가능한 우수한 내화학성
- 용기 바닥에 홈이 있어 전용 핸드렌치로 쉽게 개폐 가능 (LSMV-005 제외)
- 5ml마다 눈금이 있음(LSMV-001, LSMV-003, LSMV-005 제외)
- LK Labkorea

### Mini Vial

Cat. No.	Model	용량	Size (OD×h)	EA Price ₩
V04-1995-045	LSMV-001	3ml	∅21.8×h24 mm	13,300
V04-1995-048	LSMV-003	5ml	∅21.8×h28 mm	14,600
V04-1995-051	LSMV-005	7ml	∅21.8×h36 mm	15,700
V04-1995-054	LSMV-007	7ml	∅21.8×h36 mm	15,700
V04-1995-057	LSMV-009	12ml	∅21.8×h57.5 mm	20,400
V04-1995-060	LSMV-011	25ml	∅21.8×h110.5 mm	28,200

### Screw Cap

Cat. No.	Model	적용 Vial	EA Price ₩
V04-1995-065	LSMVCP-001	LSMV-001, LSMV-003, LSMV-005, LSMV-007, LSMV-009, LSMV-011	13,000

### Hand Wrench

Cat. No.	Model	적용 Vial	EA Price ₩
V04-1995-070	LHW-011	LSMV-001, LSMV-003, LSMV-007, LSMV-009, LSMV-011	48,200



## High Purity Teflon Sample Vial 고순도 테프론 샘플 바이알

- 재질 : PFA (Perfluoroalkoxy), 250°C 내열
- 불산 사용이 가능한 우수한 내화학성
- 샘플 저장 용도 (가열 용기로 사용 불가)
- LK Labkorea

### Sample Vials

Cat. No.	Model	용량	Size (OD×h)	EA Price ₩
V04-1995-010	LSV-001	15 ml	∅29×h39 mm	27,700
V04-1995-013	LSV-003	30ml	∅29×h70 mm	36,100
V04-1995-016	LSV-005	50ml	∅29×h117 mm	60,200

### Screw Cap

Cat. No.	Model	최대 허용 압력	적용 Vial	EA Price ₩
V04-1995-020	LSVCP-001	2 bar	LSV-001, LSV-003, LSV-005	18,300

### Stopper lid

Cat. No.	Model	적용 Vial	EA Price ₩
V04-1995-025	LSL-001	For ∅31mm Vials, LSV-001, LSV-003, LSV-005	120,400

### Dropping bottle

Cat. No.	Model	용량	Size (OD×h)	EA Price ₩
V04-1995-030	LDB-001	30ml	∅27.3×h101.5 mm	80,600

### Hand Wrench

Cat. No.	Model	Description	EA Price ₩
V04-1995-035	LHW-001	For ∅31mm Vials	48,200

