

Serum Vial / Bottle

세럼 바이알 [병]

- 재질 : Borosilicate Glass 3.3
- 캡 별도 구매
- LK Labkorea

투명 세럼 바이알 / 병

Cat. No.	Model	용량	Size	Seal 규격	PK Price ₩
V05-215-404	13-460	2 ml	Ø15×h33 mm	13 mm	100ea 31,000
V05-215-408	13-461	3 ml	Ø16.4×h36 mm		100ea 32,000
V05-215-412	13-462	5 ml	Ø22×h40 mm	20 mm	100ea 40,000
V05-215-416	13-463	10 ml	Ø24×h56 mm		100ea 47,000
V05-215-420	13-464	20 ml	Ø30×h60 mm		100ea 81,300
V05-215-424	13-465	30 ml	Ø33.5×h63 mm		100ea 70,000
V05-215-428	13-466	40 ml	Ø35×h80 mm		100ea 72,000
V05-215-432	13-467	50 ml	Ø35×h90 mm		100ea 75,000

갈색 세럼 바이알 / 병

Cat. No.	Model	용량	Size	Seal 규격	PK Price ₩
V05-215-812	12-202A	5 ml	Ø22×h40 mm	20 mm	100ea 47,000
V05-215-816	12-203A	10 ml	Ø22×h52 mm		100ea 59,000
V05-215-820	12-204A	20 ml	Ø30×h60 mm		100ea 107,000
V05-215-828	12-206A	50 ml	Ø35×h90 mm		100ea 113,000


투명 세럼 바이알 / 병, Kimble®

Cat. No.	Model	용량	Size	Seal 규격	PK Price ₩
V05-215-406	61000G-5	5 ml	Ø23×h47 mm	20 mm	288ea 350,100
V05-215-410	61000G-10	10 ml	Ø25×h52 mm		288ea 481,000
V05-215-414	61000G-20	20 ml	Ø33×h59 mm		288ea 586,000
V05-215-418	61000G-30	30 ml	Ø37×h66 mm		288ea 735,100
V05-215-422	61000G-50	50 ml	Ø43×h75 mm		288ea 904,800
V05-215-426	61000G-60	60 ml	Ø40×h90 mm		144ea 517,100
V05-215-435	61000G-100	100 ml	Ø52×h94 mm		48ea 247,900
V05-215-437	61000G-125	125 ml	Ø54×h106 mm		24ea 227,800
V05-215-439	61000G-200	200 ml	Ø65×h114 mm		36ea 204,400

- DWK Lifescience / Kimble®

대용량 세럼병, 250 ml & 500 ml

Cat. No.	Model	용량	Size	Seal 규격	CS Price ₩
V05-215-952	13-471	250 ml	Ø60×158 mm	30 mm	35ea 508,700
V05-215-954	13-472	500 ml	Ø75×190 mm		24ea 425,040



Aluminum Seal with Septa

알루미늄 씰과 셉타 일체형

Cat. No.	Model	Seal	Septa 재질	PK Price ₩	CS Price ₩
V05-214-704	5140-08	8 mm	PTFE / Red Rubber	100ea 18,900	10pk 175,800
V05-214-708	5150-08		PTFE / Silicon	100ea 35,500	10pk 331,300
V05-214-712	5160-08		PTFE / Silicon / PTFE	100ea 43,000	10pk 401,800
V05-214-716	5110-08	11 mm	10 mil PTFE	100ea 20,900	10pk 194,500
V01-214-603	5140-11		PTFE / Red Rubber	100ea 19,300	10pk 179,600
V01-214-606	5150-11		PTFE / Silicon	100ea 39,700	10pk 371,000
V01-214-609	5160-11	13 mm	PTFE / Silicone / PTFE	100ea 45,200	10pk 421,600
V01-214-684	5140-13		PTFE / Red Rubber	100ea 45,300	10pk 422,700
V01-214-687	5150-13		PTFE / Silicon	100ea 52,900	10pk 494,300
V01-214-690	5160-13	20 mm	PTFE / Silicon / PTFE	100ea 60,200	10pk 562,700
V05-214-728	5140-20		PTFE / Gray Butyl ①	100ea 122,800	10pk 1,147,700
V05-214-730	5142-20		PTFE / Gray Butyl Flange	100ea 168,200	10pk 1,571,900
V05-214-732	5150-20		PTFE / Silicon ②	100ea 77,400	10pk 722,900
V05-214-734	5150UL-20		PTFE / Silicon	100ea 104,100	10pk 972,400

- Ultra Low Bleed Septa : 5150UL-20
 - Very Sensitive Application & Multiple Injection용
 - Eliminate Particulates with no Interfering Peaks
- 추가적인 Aluminum Seal과 Septa는 Vial란의 Autosample Vial 참조
- One-piece Aluminum Seal 색상 : Silver
- JGF (J. G. Finneran)

