



Vial

Stopper for Serum and Crimp-top / Septa for Aluminum seal / Handy Plier Decapper

Stopper for Serum and Crimp-top Vial 스토퍼

• JGF (J. G. Finneran)

Stopper for 13 mm Crimp Top

Cat. No.	Model	재질	Type	PK Price ₩	CS Price ₩
V06-214-816	5001S-13G	Gray Butyl ①	Snap-On	100ea 21,000	10pk 199,500

Stopper for 20 mm


Cat. No.	Model	재질	Type	PK Price ₩	CS Price ₩
V05-214-800	5001SW-20G	Gray Butyl ②	Straight Plug	100ea 30,600	10pk 290,700
V05-214-804	5001SW-20GT	PTFE Faced Gray Butyl ③		100ea 125,500	-
V05-214-806	5001SW-20R	Red Rubber ④		100ea 42,400	10pk 395,700
V05-214-808	5001SW-20S	Clear Silicon ⑤		100ea 85,250	10pk 738,100
V05-214-810	5001SW-20V	Clear Viton®		100ea 233,640	-
V05-214-816	5001-20ST	Gray Butyl, Siliconized ⑥	Snap-On	100ea 43,300	10pk 404,700
V05-214-818	5001S-20RS	Red Rubber, Siliconized		100ea 55,600	10pk 556,000

Stopper for 30 mm Crimp Top

Cat. No.	Model	재질	Type	PK Price ₩	CS Price ₩
V06-214-830	5001SW-30G	Gray Butyl ⑦	Straight Plug	100ea 103,900	10pk 970,800
V06-214-832	5001-30ST	Gray Butyl	Snap-On	100ea 115,700	10pk 1,105,700

Lyophilization Stopper, 동결전조용 스톱퍼

Cat. No.	Model	Size	재질	Type	PK Price ₩	CS Price ₩
V05-214-814	5001L2-13	13 mm	Gray Butyl	2-Leg	100ea 22,300	10pk 207,900
V05-214-824	5001L2-20	20 mm		2-Leg	100ea 29,400	10pk 279,300
V05-214-828	5001L3-20			3-Leg ⑧	100ea 26,400	10pk 250,800



Septa for Aluminum Seal 알루미늄 씰용 셉타

Cat. No.	Model	사용 Seal	Septa 재질	PK Price ₩	CS Price ₩
V01-214-740	606550T-08	8 mm	PTFE / Silicon	100ea 21,200	10pk 197,700
V06-214-248	607550T-13	13 mm	PTFE / Silicon	100ea 38,200	10pk 356,700
V05-214-792	610050-20S	20 mm	PTFE / Silicon	100ea 64,700	10pk 604,600
V05-214-788	610040G-20		PTFE / Gray Butyl	100ea 100,000	10pk 934,600

• JGF (J. G. Finneran)



FEP/Silicon Septa for Aluminum Seal 알루미늄 씰용 FEP/Silicon 셉타

Cat. No.	Model	Thickness		GPI Finish	PK Price ₩	CS Price ₩
		FEP	Silicon Rubber			
V05-137-606	N73831A-20	0.076 mm	3.048 mm	20 or 20A	100ea 234,000	10pk 2,186,500

- Temp. Range : -60 ~ 200°C
- Headspace Autosampler에 사용되는 샘플 바이알 Sealing에 사용
- DWK Lifescience / Kimble®



Handy Plier Decapper 집게형 디캡퍼

Cat. No.	Model	사용 Seal 규격	EA Price ₩
V06-214-978	9310-08	8 mm	126,600
V06-214-980	9310-11	11 mm	126,600
V06-214-984	9310-13	13 mm	126,600
V06-214-988	9310-20	20 mm	147,500
V06-214-992	9310-30	30 mm	151,800

• JGF (J. G. Finneran)

