

Vial Rack, PP 바이알 랙

Cat. No.	Model	홀 (Ø)	홀 수	Size (w×d×h)	EA Price ₩
V05-55-102	RWV001	9.5 mm	96	125×185×22 mm	148,100
V05-55-110	RWV002	12.5 mm	50	100×190×22 mm	88,900
V05-55-118	RWV004	15.5 mm	48	94×266×28 mm	100,700
V05-55-126	RWV008	17.1 mm	90	170×330×30 mm	118,400
V05-55-134	RWV016	23.1 mm	36	91×322×28 mm	88,200
V05-55-142	RWV020	28.1 mm	50	170×330×30 mm	149,300
V05-55-150	RWV050	30.0 mm	50	170×337×30 mm	156,900

- 재질 : Polypropylene (내열 온도 : -10 ~ +140°C), Autoclavable
- 단단식으로 쌓아 올릴 수 있음 • LabPlasti®



Model	Type	홀 (Ø)	홀 수	배열	적용 바이알
RWV001	985750	9.5 mm	96	8 × 12	od8 mm Chromatography Vial
RWV002	985800	12.5 mm	50	5 × 10	od12 mm (1.8 / 2.0 ml) Autosample Vial
RWV004	868804	15.5 mm	48	4 × 12	4 ml Sample Vial
RWV008	868810	17.1 mm	90	6 × 15	8 ml Sample Vial
RWV016	868805	23.1 mm	36	3 × 12	16 ml Sample Vial, 6/10/20 ml Headspace Vial
RWV020	868806	28.1 mm	50	5 × 10	20 / 25 / 40 ml Sample Vial, 20 ml Scintillation Vial
RWV050	868808	30.0 mm	50	5 × 10	od30 mm까지의 Vial과 Bottle

Vial Box / Vial Rack, PS 바이알 보관 박스 / 바이알 랙

Cat. No.	Model	홀 (Ø)	홀 수	Type	EA Price ₩
V05-215-901	703-01	13 mm	78 (6×13)	A	25,800
V05-215-905	703-02	16 mm	60 (5×12)		25,800
V05-215-909	703-03	19 mm	50 (5×10)		25,800
V05-215-913	703-04	22 mm	32 (4×8)		25,800
V05-215-917	703-05	25 mm		B	24,000
V05-215-921	703-06	28 mm	18 (3×6)		26,400
V05-215-925	703-07	31 mm			26,400
V05-215-929	703-08	36 mm	26,400		
V05-215-933	703-11	Lid Only, w138×d258×h11.8 mm			10,500


- 재질 : PS, 내열 온도 : -30 ~ 60°C
- 구성 : Box + Holder • 뚜껑은 별도 구매
- A형 Size : w133×d253×h58.7 mm, Holder 높이 25 mm (Hold는 분리형)
- B형 Size : w133×d253×h88.7 mm, Holder 높이 35 mm (Hold는 분리형)
- Universal



Vial Storage Box 바이알 보관 박스

Cat. No.	Model	홀 수	배열	홀 내경	EA Price ₩
V05-215-941	MVC-13	84	7×12	13 mm	24,300
V05-215-945	MVC-16	60	6×10	16 mm	23,600

- 재질 : AS
- Size : w170.6×d104.4×h65 mm
- Universal



Chem-Bin® Cabinet for Vial and Closure 바이알 및 뚜껑 보관 캐비닛

Cat. No.	Model	서랍 수	Size (w×d×h)	EA Price ₩
V06-132-R68	5730-0009	9ea	432×279×279 mm	334,700

- Powder Coated Finished
- Vial 및 Closure의 보관 및 정리용
- Chemglass

