

Anaerobic Media Bottle

혐기성 미디어 바틀

Cat. No.	Model	용량	Size (Ø × h)	EA Price ₩
B08-132-110	CLS-4217-01	50 ml	46 × 88 mm	179,000
B08-132-115	CLS-4217-02	100 ml	56 × 100 mm	182,400
B08-132-120	CLS-4217-03	250 ml	70 × 138 mm	191,600
B08-132-125	CLS-4217-04	500 ml	85 × 176 mm	194,600
B08-132-130	CLS-4217-05	1 L	100 × 215 mm	210,200

- Borosilicate Glass 3.3, ASTM E438, Type I, Class A
- 메탄 생성 박테리아, Cell Culture, 신약 개발, R&D에 사용하기 적합
- Rubber Stopper and Aluminum Seal 포함, Autoclavable
- Chemglass

Optional Accessory

Cat. No.	Model	Description	PK Price ₩
T16-132-017	CLS-4209-14	Blue Butyl Rubber Stopper	100ea 154,800
V06-214-164	5100-20	Open-top Aluminum Seal	100ea 10,200

- Size : Ø20mm용



Anaerobic Bottle

혐기성 바틀

Cat. No.	Model	용량	Neck	EA Price ₩
B08-132-135	CLS-4215-03	50 ml	20 mm 1-Neck	254,800
B08-132-140	CLS-4215-05	1L		324,600
B08-132-145	CLS-4215-01	2L	20 mm 2-Neck	393,400

- Aluminum Seal 및 Stopper는 별도 구매해야 함 (1253 ~ 1255 Page 참조)
- Chemglass



Heavy-wall Hybridization Bottle

헤비월 하이브리드 병

- Heavy Wall Type으로 하이브리디제이션 인큐베이터 또는 오븐에서 사용 가능
- 직선형으로 멤브레인 배치 및 제거, Solution Decanting이 쉬움
- Red Cap : PTFE faced liner & up to 180°C, GL45
- Blue Cap : Silicone sealing ring & up to 140°C, GL45
- Dual Chamber Bottle은 양쪽에 Cap이 있으며 중앙의 유리벽이 동일한 칸으로 나눔
- Chemglass

with Red Cap

Cat. No.	Model	Size (Ø × h)	Chamber	EA Price ₩
T16-132-431	CG-1140-01	35 × 75 mm	Single	174,600
T16-132-437	CG-1140-02	35 × 100 mm		175,100
T16-132-443	CG-1140-03	35 × 150 mm		175,400
T16-132-449	CG-1140-04	35 × 250 mm		161,400
T16-132-455	CG-1140-05	35 × 300 mm		182,500
T16-132-601	CG-1140-06	35 × 150 mm	Dual	210,500
T16-132-607	CG-1140-07	35 × 300 mm		222,000

with Blue Cap

Cat. No.	Model	Size (Ø × h)	Chamber	EA Price ₩
T16-132-613	CG-1140-10	35 × 75 mm	Single	173,600
T16-132-619	CG-1140-11	35 × 100 mm		174,100
T16-132-625	CG-1140-12	35 × 150 mm		174,500
T16-132-631	CG-1140-13	35 × 250 mm		177,100
T16-132-637	CG-1140-14	35 × 300 mm		181,600
T16-132-643	CG-1140-16	35 × 150 mm	Dual	208,000
T16-132-649	CG-1140-17	35 × 300 mm		219,800



with Red Cap



with Blue Cap