

Glassware Washer with Drying System & Detergent Compartment 초차 세척기, 146 L

Cat. No.	Model	Description	EA Price ₩
W07-78-040	GW4190	Glassware Washer 146L w/ Drying System	

- Designed to work with either one to two independent and flexible washing levels
- Tank Size : w525×d490×h570 mm (146L)
- Washing up to 95°C – antibacterial thermal disinfection
- Microprocessor programmable electronic control
- Serial port RS232 for connection to external printer or PC
- Forced hot air drying system with air filter 98% DOP
- Alkaline powder detergent dispenser
- Washing chamber and internal door in AISI 316L stainless steel
- Integrated compartment for the detergents
- Detergent cabinet made of AISI 304 stainless steel 4 x 5 liters tanks
- 2 independent washing levels through telescopic guides
- 다양한 Optional Rack을 이용 실험실에서 사용되는 거의 모든 초차를 세척 가능
- Electronic door lock system for user safety
- Without trolleys
- smeg



LPV40DS



LM20DS



LB32DS



PD70I

Optional Rack

Cat. No.	Model	Description	EA Price ₩
W07-78-701	LPV40DS	Rack for 100 Pipettes	
W07-78-705	LM20DS	Universal Rack for Narrow Neck glassware 24 Position	
W07-78-709	LB32DS	Rack for Bottle	
W07-78-713	PD70I	Rack for Petridish	

Cat. No.	W07-78-040
Model	GW4190
Electronic Control	Microprocessor
Standard Program	9
Customisable Program	6
Tank Internal Temp.	From Room Temp. to 95°C
Safety Lock	Safety Thermostats, Door Interlock
Communication	RS232 Serial Port for PC or Printer
Inlet	Cold Water (Max. 6°C), Demineralized Water (<30us/cm)
Recirculation Pump	400 L/min
Drying System	Fan, 0.8 kW (Heating Element), HEPA Filter (Option)
Heating	Electrical 2.3 kW Max.
Material	Stainlss Steel (Tank : 316L, External : 304)
Outside Dimension	w900×d600×h850mm
Weight	105kg
Power	230V, 60Hz

- 더 다양한 초차 세척기는 www.labmart.co.kr 참조

