

Watch Glass / Beaker Cover 시계 접시, PP & PTFE

• 용도 : 시계 접시 또는 Beaker Cover로 사용

Sodalime Glass 재질

Cat. No.	Model	Diameter	EA Price ₩	PK Price ₩
B05-53-616	WG-G040	Ø40 mm	2,000	10ea 17,800
B05-53-624	WG-G050	Ø50 mm	2,000	10ea 17,800
B05-53-632	WG-G060	Ø60 mm	2,100	10ea 18,400
B05-53-640	WG-G070	Ø70 mm	2,200	10ea 19,300
B05-53-648	WG-G080	Ø80 mm	2,200	10ea 20,000
B05-53-656	WG-G090	Ø90 mm	2,400	10ea 21,400
B05-53-664	WG-G100	Ø100 mm	2,700	10ea 24,500
B05-53-672	WG-G120	Ø120 mm	4,200	10ea 37,600
B05-53-680	WG-G125	Ø125 mm	4,400	10ea 39,800
B05-53-688	WG-G150	Ø150 mm	4,900	10ea 44,300
B05-53-696	WG-G200	Ø200 mm	8,100	10ea 73,700
B05-53-704	WG-G250	Ø250 mm	12,600	5ea 56,900

• 재질 : Sodalime Glass, Max.250°C 내열

• Techno Glass

PP 재질

Cat. No.	Model	Material	Diameter	EA Price ₩
B05-339-001	46080-0060	PP	Ø60 mm	1,080
B05-339-003	46080-0080		Ø80 mm	1,450
B05-339-005	46080-0100		Ø100 mm	1,930
B05-339-007	46080-0120		Ø120 mm	2,410

• 재질 : PP, -10 ~ 135°C 내열

• LabPlasti®

PTFE 테프론 재질

Cat. No.	Model	Material	Diameter	EA Price ₩
B05-450-586	013.030	PTFE	Ø30 mm	5,800
B05-450-590	013.040		Ø40 mm	7,100
B05-450-594	013.050		Ø50 mm	7,600
B05-450-600	013.065		Ø65 mm	8,400
B05-450-606	013.075		Ø75 mm	9,100
B05-450-612	013.080		Ø80 mm	9,400
B05-450-618	013.100		Ø100 mm	13,000
B05-450-624	013.125		Ø125 mm	20,200
B05-450-630	013.150		Ø150 mm	31,600
B05-450-636	013.200		Ø200 mm	54,300

• 재질 : PTFE, -200°C ~ 280°C 내열

• Cowie

FEP 테프론 재질

Cat. No.	Model	Material	Diameter	EA Price ₩
B05-64-666	F1180-03	FEP	Ø50 mm	37,300
B05-64-672	F1180-05		Ø75 mm	59,400
B05-64-678	F1180-06		Ø90 mm	62,800
B05-64-681	F1180-07		Ø100 mm	78,100
B05-64-684	F1180-08		Ø125 mm	204,600


• 재질 : FEP, -200°C ~ 205°C 내열

• Flon Industry

Quartz 재질

Cat. No.	Model	Material	Diameter	EA Price ₩
B05-60-500	QH-60	Quartz	Ø60 mm	84,600
B05-60-504	QH-90		Ø90 mm	110,600
B05-60-508	QH-100		Ø100 mm	119,700
B05-60-512	QH-120		Ø120 mm	130,000
B05-60-516	QH-150		Ø150 mm	208,000

• LK Labkorea



Sodalime Glass




PP



PTFE



FEP



Quartz