

RO & UP Water Purification 순수 / 초순수 제조장치, NEW Human Series

- 실험실 및 연구소등에 널리 이용 되어지는 초순수 제조 장치임
- Water Sensor는 자동으로 Calibration되며, 적당한 수질을 유지 시키기 위해 Membrane 자동 세척기능이 있음 (시간 설정 가능)
- RS232C Port 내장 • **Certified CE, ISO9001, ISO14001**
- Human

RO & UP System, RO(순수), UP(초순수) 겸용 사용

배송비 별도

Cat. No.	Model	Capacity	Type	EA Price ₩
W02-78-311	New Human Power I S	15 L/hr	Scholar	
W02-78-318	New Human Power I UV		Schoalr-UV	
W02-78-325	New Human Power I B		Bio	
W02-78-332	New Human Power I IN		Integrate	
W02-78-339	New Human Power II S	25 L/hr	Scholar	
W02-78-346	New Human Power II UV		Schoalr-UV	
W02-78-353	New Human Power II B		Bio	
W02-78-360	New Human Power II IN		Integrate	
W02-78-367	New Human Power III S	35 L/hr	Scholar	
W02-78-374	New Human Power III UV		Schoalr-UV	
W02-78-381	New Human Power III B		Bio	
W02-78-388	New Human Power III IN		Integrate	

• Size : w400×d567×h510 mm

RO System, RO(Reverse Osmosis) 전자동 순수 제조 장치

배송비 별도

Cat. No.	Model	Capacity	Type	EA Price ₩
W02-78-450	New Human RO180	15 L/hr	-	
W02-78-457	New Human RO280	25 L/hr		
W02-78-464	New Human RO380	35 L/hr		

• Size : w348×d560×h518 mm

UP System, UP(초순수) 제조 장치

배송비 별도

Cat. No.	Model	Capacity	Type	EA Price ₩
W02-78-507	New Human UP900 S	1.5 L/min	Scholar	
W02-78-514	New Human UP900 UV		Scholar-UV	
W02-78-521	New Human UP900 B		Bio	
W02-78-528	New Human UP900 IN		Integrate	

• Size : w348×d560×h518 mm

공통 Specification

Model		New Human Power -				New Human RO-				New Human UP900-			
Water Quality Display	Feed Water	0 to 990 us/cm				-				0 to 250 us/cm			
	RO	0.2 to 250.0 us/cm				-				-			
	UP	Up to 18.3 MΩ-cm				-				Up to 18.3 MΩ-cm			
	Temp.	0 to 99.9°C											
Controller		Microprocess One Touch Technology											
Display		16 X 2 Characters Backlight LCD Digital Display											
Power		230V, 50/60 Hz, 60W											
Type	Scholar	Scholar-UV	Bio	Intergrate	-	Scholar	Scholar-UV	Bio	Intergrate				
0.2um Final Filter	O	O	O	O		O	O	O	O				
258/185 UV Lamp	X	O	X	O		X	O	X	O				
MW 5000 UF Filter	X	X	O	O		X	X	O	O				
Dispenser	-	-	-	-		-	-	-	-				

Spare Accessories for New Human Series

Cat. No.	Model	교체주기	Type	EA Price ₩
W02-78-904	HMC-10PF	약 1개월	Pre Filter	3,200
W02-78-030	New Human RO Pack I	약 18개월	New Human Power I & RO180	291,000
W02-78-034	New Human RO Pack II	약 18개월	New Human Power II & RO 280	396,000
W02-78-038	New Human RO Pack III	약 18개월	New Human Power III & RO 380	465,000
W02-78-042	New Human UP Pack S	약 18개월	Scholar	429,000
W02-78-044	New Human UP Pack B	약 18개월	Bio	429,000
W02-78-046	New Human UP Pack S-UV	약 18개월	Scholar-UV	429,000
W02-78-048	New Human UP Pack I	약 18개월	Intergrate	429,000
W02-78-920	HMC-0.2CF	약 8개월	0.2um Capsule Filter	71,000
W02-78-930	HMC-254UV	약 2000시간	Scholar-UV & Intergrate	332,000
W02-78-941	HMC-MW5000L	약 8개월	Bio & Intergrate	758,000



Capsule Filter

RO & UP System



CE

RO System



CE

UP System



CE