

I Glass Narrow Neck Bottle 세구병, with Teflon Lined Cap

- 1개의 Cell에 1개의 Bottle이 들어있는 매우 편리한 청결 포장
- Green Thermoset F217 Teflon Lined Cap일체형
- ! Do not use with oven and autoclave
- Qorpak

Clear Narrow Neck Bottle

Cat. No.	Model	용량	Cap	Size	PK Price ₩
B08-134-005	GLC01082	30 ml	20-400	Ø31×h79 mm	48ea 128,800
B08-134-010	GLC01099	60 ml		Ø39×h94 mm	24ea 77,800
B08-134-015	GLC01123	125 ml	22-400	Ø48×h112 mm	24ea 87,100
B08-134-020	GLC01151	250 ml	24-400	Ø60×h137 mm	24ea 111,900
B08-134-025	GLC01185	500 ml	28-400	Ø75×h168 mm	12ea 69,000
B08-134-030	GLC01217	1 L	33-400	Ø94×h210 mm	12ea 98,300

Amber Narrow Neck Bottle

Cat. No.	Model	용량	Cap	Size	PK Price ₩
B08-134-095	GLC01904	60 ml	20-400	Ø39×h94 mm	24ea 79,000
B08-134-100	GLC01926	125 ml	22-400	Ø48×h112 mm	24ea 88,900
B08-134-105	GLC01947	250 ml	24-400	Ø60×h138 mm	24ea 105,500
B08-134-110	GLC01977	500 ml	28-400	Ø75×h168 mm	12ea 77,800
B08-134-115	GLC02004	1 L	33-400	Ø94×h210 mm	12ea 102,300



I Glass Narrow Neck Bottle 세구병, with PE Cone Lined Cap

- 1개의 Cell에 1개의 Bottle이 들어있는 매우 편리한 청결 포장
- Black Phenolic Cap with PE Cone
- ! Do not use with oven and autoclave
- JGF (J. G. Finneran)

Clear Narrow Neck Bottle

Cat. No.	Model	용량	Cap	Size	PK Price ₩
B08-214-541	D0149-1/2	15 ml	18-400	Ø25×h68 mm	24ea 59,900
B08-214-543	D0149-1	30 ml		Ø31×h79 mm	24ea 58,300
B08-214-545	D0149-2	60 ml	20-400	Ø39×h94 mm	24ea 65,100
B08-214-547	D0149-4	125 ml	22-400	Ø48×h112 mm	24ea 68,300
B08-214-549	D0149-8	250 ml	24-400	Ø60×h136 mm	12ea 36,400
B08-214-551	D0149-16	500 ml	28-400	Ø75×h168 mm	12ea 58,600
B08-214-553	D0149-32	1 L	33-400	Ø94×h206 mm	12ea 69,900

Amber Narrow Neck Bottle

Cat. No.	Model	용량	Cap	Size	PK Price ₩
B08-214-581	D0161-1/2	15 ml	18-400	Ø25×h68 mm	24ea 59,900
B08-214-583	D0161-1	30 ml		Ø31×h79 mm	24ea 58,300
B08-214-585	D0161-2	60 ml	20-400	Ø39×h94 mm	24ea 65,400
B08-214-587	D0161-4	125 ml	22-400	Ø48×h112 mm	24ea 71,800
B08-214-589	D0161-8	250 ml	24-400	Ø60×h136 mm	12ea 41,700
B08-214-591	D0161-16	500 ml	28-400	Ø75×h168 mm	12ea 58,600
B08-214-593	D0161-32	1 L	33-400	Ø94×h206 mm	12ea 76,700

