

Clear Tall Straight Side Round Bottle / Jar
장형 대 광구병, with F217 Foam Lined

Cat. No.	Model	용량	Cap	Size	PK Price ₩
B08-214-511	D0098-1	30 ml	33-400	Ø34×h68 mm	24ea 59,200
B08-214-513	D0098-2	60 ml	38-400	Ø42×h83 mm	24ea 67,500
B08-214-515	D0098-4	125 ml	48-400	Ø51×h102 mm	24ea 64,600
B08-214-517	D0098-8	250 ml	58-400	Ø62×h127 mm	12ea 37,200
B08-214-519	D0098-16	500 ml	70-400	Ø76×h145 mm	12ea 53,600

- 1개의 Cell에 1개의 Bottle이 들어있는 매우 편리한 청결 포장
- White PP Cap with F217 Foam Lined
- ! Do not use with oven and autoclave
- JGF (J. G. Finneran)



Clear Square Tablet Bottle / Jar
사각 태블릿 광구병

Cat. No.	Model	용량	Cap	Size	PK Price ₩
B08-134-470	GLC01275	30 ml	33-400	Ø34×h60 mm	48ea 163,000
B08-134-475	GLC01302	60 ml		Ø42×h85 mm	24ea 93,900

- 1개의 Cell에 1개의 Bottle이 들어있는 매우 편리한 청결 포장
- Green Thermoset F217 Teflon Lined Cap일체형
- ! Do not use with oven and autoclave
- Qorpak



Wide Mouth French Square Bottle / Jar
사각 대 광구병, w/ Teflon Lined Cap

Cat. No.	Model	용량	Cap	Size	PK Price ₩
B08-134-200	GLC01241	15 ml	20-400	Ø24×h62 mm	48ea 151,600
B08-134-204	GLC01260	30 ml	24-400	Ø30×h72 mm	48ea 173,900
B08-134-208	GLC01284	60 ml	28-400	Ø37×h86 mm	48ea 178,000
B08-134-212	GLC01311	125 ml	33-400	Ø45×h112 mm	24ea 114,700
B08-134-216	GLC01336	250 ml	43-400	Ø51×h137 mm	24ea 159,900
B08-134-220	GLC01364	500 ml	48-400	Ø62×h167 mm	24ea 189,200
B08-134-222	GLC01388	1 L	58-400	Ø83×h202 mm	12ea 189,600

- 1개의 Cell에 1개의 Bottle이 들어있는 매우 편리한 청결 포장
- Green Thermoset F217 Teflon Lined Cap일체형
- ! Do not use with oven and autoclave
- Qorpak



Glass Jug With Teflon Lined Cap
글라스 저그

- Qorpak

With Green Thermoset F217 Teflon Lined Cap

Cat. No.	Model	색상	용량	Cap	PK Price ₩
B08-134-502	GLC01408	Clear ①	2 L	38-400	6ea 114,000
B08-134-509	GLC01429		3.8 L		4ea 100,800
B08-134-511	GLC02043	Amber ②			6ea 135,300

- Green Thermoset F217 Teflon Lined Caps일체형
- ! Do not use with oven and autoclave

With Black PP Cap & Teflon Disc

Cat. No.	Model	색상	용량	Cap	PK Price ₩
B08-134-512	GLC02062	Amber ③	2.5 L	38-430	6ea 119,900
B08-134-518	259531		4 L	38-439	6ea 173,900

- ! Do not use with oven and autoclave

