

I Screw Cap for Bottle


스크류 캡

• Qorpak

Teflon (PTFE) – Lined Green Thermoset Screw Cap

Cat. No.	Model	Cap	PK Price ₩		CS Price ₩	
B08-134-C01	C-00544	8-425	100ea	79,300	5pk	370,500
B08-134-C02	C-00546	13-425	100ea	58,700	5pk	274,300
B08-134-C03	C-00548	15-425	100ea	66,900	5pk	312,400
B08-134-C04	C-00550	18-400	100ea	79,400	5pk	397,100
B08-134-C05	C-00552	20-400	100ea	97,300	5pk	454,500
B08-134-C06	C-00554	22-400	100ea	126,500	5pk	591,000
B08-134-C07	C-00556	24-400	100ea	131,000	5pk	654,900
B08-134-C08	C-00558	28-400	100ea	156,900	5pk	733,100
B08-134-C09	C-00560	33-400	100ea	217,400	5pk	1,015,900
B08-134-C10	C-00562	38-400	100ea	288,200	5pk	1,346,600
B08-134-C11	C-00066	43-400	100ea	341,100	5pk	1,593,700
B08-134-C12	C-00068	45-400	100ea	369,800	5pk	1,728,000
B08-134-C13	C-00071	48-400	100ea	419,100	5pk	1,958,200
B08-134-C14	C-00073	53-400	100ea	497,700	5pk	2,325,300
B08-134-C15	C-00075	58-400	100ea	564,900	5pk	2,639,400
B08-134-C16	C-00077	63-400	72ea	515,800	4pk	1,928,000
B08-134-C17	C-00079	70-400	72ea	608,200	2pk	1,136,700
B08-134-C18	C-00061	83-400	72ea	956,500	2pk	1,787,800
B08-134-C19	C-00063	89-400	72ea	1,034,500	2pk	1,933,500

• 오븐 건조 및 Autoclave 불가



Poly-Seal Cone Lined Black Phenolic Cap

Cat. No.	Model	Cap	PK Price ₩		CS Price ₩	
V05-134-C32	C-00111	13-425	144ea	128,400	4pk	479,800
V05-134-C35	C-00125	15-425	144ea	132,400	4pk	494,900
V05-134-C37	C-00135	18-400	144ea	134,500	4pk	502,500
B08-134-C31	C-00148	20-400	144ea	136,800	4pk	511,100
B08-134-C32	C-00161	22-400	144ea	140,600	4pk	525,500
B08-134-C33	C-00176	24-400	144ea	142,600	4pk	532,800
B08-134-C34	C-00195	28-400	144ea	151,600	4pk	566,500
B08-134-C35	C-00212	33-400	144ea	187,000	4pk	698,800
B08-134-C36	C-00231	38-400	144ea	231,300	4pk	864,700



White Rubber Lined Black Phenolic Cap

Cat. No.	Model	Cap	PK Price ₩		CS Price ₩	
B08-134-C37	C-00475	20-400	144ea	149,700	4pk	559,600
B08-134-C38	C-00477	22-400	144ea	168,300	4pk	629,000
B08-134-C39	C-00481	24-400	144ea	166,900	4pk	667,200
B08-134-C40	C-00487	28-400	144ea	193,200	4pk	722,000
B08-134-C41	C-00489	33-400	144ea	255,500	4pk	954,900
B08-134-C42	C-00491	38-400	144ea	330,600	4pk	1,235,800
B08-134-C43	C-00493	43-400	144ea	405,200	4pk	1,514,700

• Autoclave



- Black Phenolic PTFE Faced Rubber Lined Cap은 1249 Page 참조
- Black Phenolic Cemented-In Rubber Lined Cap은 1249 Page 참조
- White Urea PTFE Faced Foam Baked Rubber Lined Cap은 1249 Page 참조

